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**Please read the instruction manual carefully before using the Spengler medical device and its accessories, and keep this manual for future reference.**

## **Spengler Otoscope**

### **Product contents**

- 1 Otoscope
- Single-use ear specula, grey, 10 pcs each: 2.5 mm and 4.0 mm
- 1 Instruction manual
- 2 Alkaline batteries (AAA/LR03)



### **Intended purpose**

Spengler otoscopes are intended to be used for examining the external auditory canal. The instrument can also be used for other lighting purposes during non-invasive examinations.

### **Contraindications**

- Avoid direct contact with mucous membranes or injured skin.
- Do not shine the light in the eyes.

## General warnings

-  When cleaning/disinfecting devices that come into contact with patients, make sure that appropriate cleaning agents/disinfectants are used.
-  Make sure to use personal protective equipment and follow all safety instructions / warnings.
-  Always remove the batteries prior to cleaning/disinfection.
-  Make sure that no liquid enters the unit.
-  Only use Spengler medical devices with Spengler replacement parts and accessories. The use of replacement parts and accessories that are not listed in the instruction manual can reduce safety.
-  During use, metal components near the instrument head can become warm. This is especially the case when the device has been turned on for a long period of time. Do not leave the device turned on if it is not being used.
-  During normal use of the device, the temperature of touchable parts can reach temperatures up to 48°C.
-  Do not service Spengler medical devices while using them on patients. Be mindful of all significant risks of reciprocal interference, which may occur when Spengler medical devices are used during special examinations or treatments.
-  To prevent and minimise electromagnetic disturbances between the Spengler medical device and other devices, please follow the EMC Guidelines for this device and the instructions of other devices.
-  It is possible for the product to switch on unintentionally when it is turned off as a result of an electrostatic malfunction. The product can be turned off at any time manually or it will turn off automatically after 3 minutes.

 If a patient has a warm sensation in the areas of contact, the diagnostic examination must be temporarily stopped and the Spengler device should be turned off to cool off.

 All serious incidents that occur in connection with the product must be reported to the manufacturer and the competent authority of the member state in which the user and/or the patient is registered.

 Do not use Spengler otoscopes to shine the light in a patient's eyes.

 Ear canal examinations may only be carried out with an ear specula attached.

 NEVER use the device without an ear specula. When in use, the ear specula may not be left inside the ear (skin contact) for longer than 1 minute.

 The pressure build-up during a pneumatic test must be done carefully (excess pressure hazard and resulting damage).

 Caution: The slit for attaching the ear specula has sharp edges, and there is risk of injury.

---

## General instructions for use

The head of the medical device can be detached from the handle by pulling it vertically.

The head can be put back on the handle in reverse order by pushing the mounting adjustment vertically onto the handle with locking pins (5) until you can feel it lock into place.

## Instructions for Spengler instrument handle

The Spengler instrument handle is only intended for use with Spengler instruments. The instruments may only be used by medically trained personnel.

### Instrument connection

Mounting adjustment with locking grooves (1).

### To turn on

Press the push-button switch (2) once.

### To adjust the brightness

To make it brighter: Hold down the push-button switch (2) when the device is in switched-on mode. When the desired level of brightness has been reached, let go.


To make it less bright: Hold down the push-button switch (2) again. When the desired level of brightness has been reached, let go.

### To turn off the device

If the medical device is on, push the switch (2) once. The medical device also has an automatic shut-off function. It will turn off the device after 3 minutes.

### To change the batteries

To change the battery, open the battery compartment lid (3) by turning it to the left. **CAUTION:** Make sure to put the poles of the new batteries in the correct position (the plus pole of the battery faces the head of the medical device and the minus pole faces the battery compartment lid). Batteries placed incorrectly can damage the Spengler medical device and lead to expiry of warranty claims.

Put the battery compartment lid (3) back on and close it by turning it to the right. Do not change the batteries when in the vicinity of patients.

### **Troubleshooting**

If the light of the device flickers during operation, please replace the batteries contained in the device.

### **Cleaning**

Clean the device with a soft, lint-free towel; if needed, moistened with a little alcohol. The Spengler medical device may not be submerged in liquid.

### **Disinfection**

Use a disinfectant wipe to disinfect or hygienically clean the external surfaces of device. Disinfectants can be used that have been approved for medical devices made of metal, plastic and glass. Follow the information provided by the manufacturer of the cleaning agent.

In accordance with RKI recommendations (Robert Koch Institute, Germany), the disinfectant must fall under the AB classification range (verifiably bactericide, fungicide and virucide) in the RKI definition list.

### **Sterilisation**

No sterilisation is permitted.

### **Disposal**

The product must be taken to a separate waste collection site for electrical and electronic equipment.

### **Additional instructions for using the otoscope**

#### Attachment of the ear specula:


Spengler otoscopes may only be placed in the ear canal with an ear specula attached.

The ear specula should be attached to the tip of the otoscope (4) with the locking pin by turning it slightly clockwise.

An ear specula does not have to be placed on the otoscope if it is being used as a light for another purpose, as long as it does not come into direct contact with the skin or mucous membranes.

### Magnifying function:

The magnifying glass in the field of view provides 3x object magnification. The magnifying glass can be turned clockwise until the end stop for instrumentation.

EN

### Use for a pneumatic otoscopy test:

Spengler otoscope can be connected to a pneumatic head with a rubber bulb and tubing for conducting a pneumatic test (6).

### Cleaning

 Ear speculas may only be used once. Treatment can lead to damage on the ear specula.

The magnifying glass must not be cleaned/disinfected with cleaning agents containing alcohol as there is a risk of clouding. Only use products that are suitable for plastic lenses.

---

### **General technical specifications**


|                                    |   |
|------------------------------------|---|
| Nominal voltage                    | 2.5 V   |
| Ampacity                           | 23mA  |
| Applied part                       | BF  |
| Protection class                   | Class II  |
| Protection rating                  | IP20  |
| Ambient conditions                 | 5° C – 40° C<br>10-90% rel. humidity<br>700-1070 hPa atmospheric pressure |
| Warehouse/<br>transport conditions | -10° and +45°C.<br>10% - 90% rel. humidity                                |
| Supply voltage                     | 2 AAA Alkaline Batteries (2.5 V variant)                                  |

## Applicable standards

IEC 60601-1  
IEC 60601-1-2  
IEC 60601-1-6

IEC 60601-1-9  
DIN EN ISO 10993-1  
IEC 62471

## Explanation of symbols

|  |   |
|--|---|
|    | Applied part type BF  |
|    | Protection class II   |
|    | Follow the instruction manual   |
|    | Stand-by (on/off button)  |
|    | Do not reuse  |
|    | Manufacturer  |
|    | Production date   |
|    | Distributor   |
|    | Caution   |
|    | The product must be taken to a separate waste collection site for electrical and electronic equipment.                    |
|    | Conforms with EU regulations  |
|  | This instrument is a medical device   |
|  | Conforms with UK regulations  |
|  | Authorised agent for medical products in Switzerland:<br>Inspira GmbH<br>Thunstrasse 64<br>CH-3110 Münsingen, Switzerland |

## Warranty conditions

We grant a warranty period of 3 years for this product starting from the date of purchase. The purchase date must be proved with a receipt. Defects due to material or production flaws will be remedied within the warranty period free of charge. The warranty period is not extended for the entire product as a result of service provided under the warranty, but rather only for the replaced parts.

Deterioration through wear and tear, transport damage and any damage caused by improper handling (e.g. not following the instruction manual) or due to interference by unauthorised persons are excluded from the warranty. The warranty does not give rise to any claims for damages against us.

The statutory claims for defects of the buyer pursuant to Section 437 of the German Civil Code (BGB) are not restricted.

In case of a claim under the warranty, the device must be sent to the manufacturer together with the original receipt:

LUXAMED GmbH & Co. KG  
Daniel-Weil-Straße 3  
89143 Blaubeuren  
GERMANY

## Customer service

Warranty and repair work must be carried out by trained and authorised personnel. The device may not be modified without the manufacturer's permission. Send the device carefully packed and with sufficient postage to your authorised dealer or directly to:

LUXAMED GmbH & Co. KG  
Abteilung Service  
Daniel-Weil-Straße 3  
89143 Blaubeuren  
GERMANY  
Tel: + 49 (0) 7344 92905-0

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# EMC Guidelines

Guidance and manufacturer declaration - electromagnetic emissions.

## Emission

| Test                       | Limit                      | Electromagnetic environment - guidance  |
|----------------------------|----------------------------|---|
| Conducted emission         | CISPR 11, Group 1, Class B | Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| Radiated emission          | CISPR 11, Group 1, Class B |   |
| Harmonic current emissions | IEC 61000-3-2, Class A     | Device is suitable for direct connection to the public low-voltage power supply network.  |

## Immunity test levels

| Test                                    | Compliance level  | Electromagnetic environment - guidance   |
|---|---|--|
| Electrostatic Discharge (IEC 61000-4-2) | Contact Discharge: $\pm 8$ kV<br>Air Discharge: $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Radiated RF EM field (IEC 61000-4-3)    | 80-2700 MHz; 1kHz AM 80 %; 10 V/m   | <p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>d = 1.2\sqrt{P}</math> for 80 MHz to 800 MHz<br/> <math>d = 2.3\sqrt{P}</math> for 800 MHz to 2,7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> |

|  |  |   |
|--|--|---|
| Proximity fields from RF wireless communications equipment (IEC 61000-4-3) | 385 MHz; Pulse Modulation: 18 Hz; 27 V/m<br>450 MHz, FM + 5 Hz deviation: 1 kHz sine; 28 V/m<br>710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m<br>810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m<br>1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m<br>2450 MHz; Pulse Modulation: 217 Hz; 28 V/m;<br>5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance 30 cm. |
|--|--|---|

|   |   |  |
|---|---|--|
| Electrical fast transients / bursts (IEC 61000-4-4) | Power lines: 2 kV; 100 kHz repetition frequency<br>Signal lines: 1 kV; 100 kHz repetition frequency | Mains power quality should be that of a typical environment. |
| Surges (IEC 61000-4-5)                              | L-N: 1kV at 0°, 90°, 180°, 270°<br>L-PE, N-PE: 2 kV at 0°, 90°, 180°, 270°                          | Mains power quality should be that of a typical environment. |

|  |  |   |
|--|--|---|
| <p>Conducted disturbances induced by RF fields (IEC 61000-4-6)</p> | <p>0.15-80 MHz; 1kHz AM 80 %; 3 Vrms, 6 Vrms in ISM and amateur radio band</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance <math>d = 1.2\sqrt{P}</math> for 150 kHz to 80MHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> |
| <p>Rated power frequency magnetic fields (IEC 61000-4-8)</p>       | <p>30 A/m, 50 Hz</p>   | <p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>  |

|  |  |  |
|--|--|--|
| <p>Voltage dips / Voltage interruptions (IEC 61000-4-11)</p> | <p>0 % UT for 0.5 cycle at 0°, 45°,90°,135°,180°,225°,270°,315°<br/> 0 % UT for 1 cycle at 0°<br/> 70 % UT for 25/30 cycles at 0°<br/> 0 % UT for 250/300 cycles at 0°</p> | <p>Mains power quality should be that of a typical environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or battery.</p> |
|--|--|--|

**[www.spengler.fr](http://www.spengler.fr)**

C1-1010\_Spengler.D.2

**UK  
CA**

**CE**