touchecg

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Bibliography

1. Pipberger HV, Arms RJ, Stallmann FW. *Automatic screening of normal and abnormal electrocardiograms by means of digital electronic computer*. Proc Soc Exp Biol Med. 1961;106:130-132.
2. Caceres CA, Steinberg CA, Abraham S, Carbery WJ, McBride JM, Tolles WE, Rikli AE. *Computer extraction of electrocardiographic parameters*. Circulation. 1962;25:356-362.
3. Macfarlane PW. *A modified axial lead system for orthogonal lead electrocardiography*. Cardiovasc Res. 1969;3:510-515.
4. Macfarlane PW. *ECG waveform identification by digital computer*. Cardiovasc Res. 1971;5:141-146.
5. Macfarlane PW, Lorimer AR, Lawrie TD. *Normal ranges of modified axial lead system electrocardiogram parameters*. Br Heart J. 1971;33:258-265.
6. Macfarlane PW, Lorimer AR, Lawrie TD. *3 and 12 lead electrocardiogram interpretation by computer. A comparison on 1093 patients*. Br Heart J. 1971;33:266-274.
7. Macfarlane PW, Cawood HT, Taylor TP, Lawrie TD. *Routine automated electrocardiogram interpretation*. Biomed Eng. 1972;7:176-180.
8. Macfarlane PW, Watts MP, Peden J, Lennox G, Lawrie TD. *Computer-assisted ECG interpretation*. Br J Clin Equip. 1976;1:61-70.
9. Taylor TP, Macfarlane PW, Lawrie TD. *Arrhythmia interpretation by digital computer*. In: Schubert E, ed. Neue Ergebnisse der Electrocardiologie II. Berlin: Humboldt University; 1974:243-245.
10. Macfarlane PW, Peden A, Podolski M, Lawrie TD. *A new 12 lead ECG diagnostic computer program*. In: Ueda H, et al, eds. Recent Advances in Electrocardiology. Jpn Heart J. 1982;23(suppl.1):667-670.
11. Macfarlane PW, Podolski M, Watts MP, Shoat D, Macfarlane DK, Lawrie TD. *The new Glasgow system*. In: Willems JL, van Bommel JH, Zyweitz C, eds. Computer ECG Analysis: Towards Standardization. Amsterdam: North Holland; 1986:31-36.
12. Macfarlane PW. *Computer interpretation of cardiac rhythm*. In: Willems JL, van Bommel JH, Zyweitz C, eds. Computer ECG Analysis: Towards Standardization. Amsterdam: North Holland; 1986:279-284.
13. Macfarlane PW, Devine B, Latif S, McLaughlin S, Shoat DB, Watts MP. *Methodology of ECG interpretation in the Glasgow program*. Meth Inf Med. 1990;29:354-361.
14. Macfarlane PW, Coleman EN, Devine B, Houston A, McLaughlin S, Aitchison TC, Pomphery EO. *A new 12-lead pediatric ECG interpretation program*. J Electrocardiol. 1990;23(suppl):76-81.
15. Yang TF, Devine B, Macfarlane PW. *Artificial neural networks for the diagnosis of atrial fibrillation*. Med Biol Eng Computing. 1994;32:615-619.
16. Yang TF, Devine B, Macfarlane PW. *Deterministic logic versus software-based artificial neural networks in the diagnosis of atrial fibrillation*. J. Electrocardiol. 1993;26(suppl):90-94.
17. McLaughlin SC, Aitchison TC, Macfarlane PW. *Improved repeatability of 12-lead ECG analysis using continuous scoring techniques*. J Electrocardiol. 1993;26(suppl):101-107.
18. McLaughlin SC, Aitchison TC, Macfarlane PW. *Methods for improving the repeatability of automated ECG analysis*. Methods Inf Med. 1995;34:272-282.
19. McLaughlin SC, Aitchison TC, Macfarlane PW. *The value of the coefficient of variation in assessing repeat variation in ECG measurements*. Eur Heart J. 1998;19:342-351.
20. Macfarlane PW, Lawrie TD, eds. *Comprehensive Electrocardiology*. New York: Pergamon Press; 1989.
21. Macfarlane PW, McLaughlin SC, Devine BD, Yang TF. *Effects of age, sex and race on ECG interval measurements*. J Electrocardiol. 1994;27(suppl):14-19.
22. Willems JL, Abreu-Lima C, Arnaud P, et al, incl. Macfarlane PW. *The diagnostic performance of computer programs for the interpretation of electrocardiograms*. New Engl J Med. 1991;325:1767-1773.
23. Morrison S, Macfarlane PW. *Computer detection of atrial flutter*. Annals of Noninvasive Electrocardiology. 2000;5:358-364.

24. Macfarlane PW. *Age and sex related normal limits of ST amplitude*. J Electrocardiol. 2001;34(suppl):235-241.
25. Macfarlane PW, Browne D, Devine B, Clark E, Miller E, Seyal J, Hampton D. *Modification of ACC/ESC criteria for acute myocardial infarction*. J Electrocardiol. 2004;37(suppl):98-103.
26. Macfarlane PW, Browne D, Devine B, Clark E, Miller E, Seyal J, Hampton D. *Effect of age and gender on diagnostic accuracy of ECG diagnosis of acute myocardial infarction*. In: A. Murray, ed. Computers in Cardiology. 2004;31:165-168.
27. Macfarlane PW, Hampton DR, Clark E, Devine B, Jayne CP. *Computer and cardiologist diagnosis of ST-elevation myocardial infarction*. J Electrocardiol. 2007;40(suppl 1):S32-S33.
28. Sgarbossa EB, Pinski SL, Barbagelata A, Underwood DA, Gates KB, Topal EJ, Califf RM, Wagner GS. *Electrocardiographic diagnosis of evolving acute myocardial infarction in the presence of left bundle-branch block*. N Engl J Med. 1996;334:481-487.
29. Tabas JA, Rodriguez RM, Seligman HK, Goldschlager NF. *Electrocardiographic criteria for detecting acute myocardial infarction in patients with left bundle branch block: a meta-analysis*. Ann Emerg Med. 2008;52:329-336.
30. Macfarlane PW, Devine B, Clark E. *The University of Glasgow (Uni-G) ECG analysis program*. Computers in Cardiology. 2005;32:451-454.
31. The CSE Working Party (incl. Macfarlane PW). *Recommendations for measurement standards in quantitative electrocardiography*. Eur Heart J. 1985;6:815-825.
32. Bailey JJ, Berson AS, Garson A, Horan LG, Macfarlane PW, Mortara D, Zywiets C. *Recommendations for standardization and specifications in automated electrocardiography: bandwidth and digital signal processing*. Circulation. 1990;81:730-739.
33. Willems JL, Arnaud P, van Bommel JH, et al, incl. Macfarlane PW. *A reference database for multilead electrocardiographic computer measurement programs*. J Am Coll Cardiol. 1987; 10:1313-1321.
34. Willems JL, Arnaud P, van Bommel JH, Degani R, Macfarlane PW, Zywiets C, for the CSE Working Party. *Common standards for quantitative electrocardiography: goals and main results*. Methods Inf Med. 1990;29:263-271.
35. Wagner G, Lim T, Gettes L, Gorgels A, Josephson M, Wellens H, Anderson S, Childers R, Clemmensen P, Kligfield P, Macfarlane P, Pahlm O, Selvester R. *Consideration of pitfalls in and omissions from the current ECG standards for diagnosis of myocardial ischemia/infarction in patients who have acute coronary syndromes*. Cardiol Clin. 2006;24:331-342.
36. Kligfield P, et al, incl. Macfarlane PW. *Guidelines for the standardization and interpretation of the electrocardiogram, part 1: the electrocardiogram and its technology*. Circulation. 2007;115:1306-1324.
37. Thygesen K, Alpert JS, White HD; Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. *Universal definition of myocardial infarction*. Circulation. 2007;116:2634-2653. Also in: J Am Col Cardiol . 2007;50:2173-2195. Also in: Eur Heart J. 2007;28:2525-2538.

Siège social

Via De Zinis, 6
38011 Cavareno (TN), Italy
T. +39 0463 850125
F. +39 0463 850088

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Via F.lli Bronzetti, 8
20129 Milano, Italy
T. +39 02 94750470
F. +39 02 94750471