
Spengler

CARDIOMATE 3[®] & CARDIOMATE 6[®]

Électrocardiographe
Electrocardiograph
Electrocardiografo
Elektrokardiograph
جهاز تخطيط القلب

CE
0123

FR	Manuel utilisateur	P.4
EN	User Manual	P.52
ES	Manual del usuario	P.100
DE	Benutzerhandbuch	P.148
AR	دليل المستخدم	P.194

WWW.SPENGLER.FR

Avis

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The 'Barcode' configuration window has a title bar with 'Barcode' and 'Operate : Select Switch'. It contains a table with three columns: 'Item', 'Start Address', and 'End Address'. The 'Item' column lists various fields: ID, First Name, Last Name, Sex, Year Of Birth, Month Of Birth, Day Of Birth, Male Code, and Female Code. The 'Start Address' and 'End Address' columns show the corresponding address ranges for each item. At the bottom, there are two buttons: 'Confirm F1' and 'Cancel F2'.

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The 'Patient' details window has a title bar with 'Patient' and 'Operate : Select Switch'. It contains a form with various fields for patient information. The 'ID' field is highlighted with a blue border and contains the value '97870401'. Below the 'ID' field, there is a barcode and the number '9 787040 195835 >'. The form includes fields for 'Admission No.', 'First Name', 'Outpatient No.', 'Last Name', 'Department', 'Age', 'Weight', 'Date of Birth', 'Height', 'Sex', 'BP (S/D)', 'Check Dr.', and 'Diagnosis Dr.'. At the bottom, there are five buttons: 'Confirm F1', 'Cancel F2', 'Clear F3', 'Copy F4', and 'Worklist F5'. An 'En' button is also present in the bottom right corner.

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
The 'Information setup' screen has a title bar with 'Information setup' on the left and 'Operate : ⬅ Select ➡ Switch' on the right. The main area is black with white text. It contains two input fields: 'Hospital name:' followed by a white rectangular box, and 'Device name:' followed by a white rectangular box. At the bottom, there are two buttons: 'Confirm F1' and 'Cancel F2'. A small 'en' button is located to the right of the 'Device name' field.

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The 'Admin configuration' screen has a title bar with 'Admin configuration' on the left and 'Operate : ⬅ Select ➡ Switch' on the right. The main area is black with white text. It contains several settings, each with radio button options:

- ID generated by: ☐ Acc. ☒ Time ☐ User input
- Upload: ☐ Auto ☒ Manual
- Upload Method: ☐ FTP ☐ Shared Directory ☒ DICOM
- File Format: ☐ FDA-XML ☒ EM-XML ☐ SCP ☐ BKG ☐ DICOM
- Picture Type: ☐ Don't save ☒ PDF ☐ JPG ☐ BMP
- Picture Direction: ☐ Horizontal ☒ Vertical
- Picture ECG Grid: ☐ Off ☐ Dot ☒ Line
- Picture Grid Color: ☐ Gray ☐ Blue ☐ Red
- Patient Query Protocol: ☐ Off ☐ HL7 ☒ Worklist
- Data Mode: ☐ ECG mode ☒ Demo ☐ Calibration

At the bottom, there are five buttons: 'Language F1', 'Logo F2', 'Update F3', 'Barcode F4', and 'More F5'.

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Fabricant : Shenzhen ECGMAC Medical Electronics Co., Ltd

Adresse de l'usine : 2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District, 518102 Shenzhen, Chine

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Site web : www.wellkang.ltd.uk, www.CE-marking.eu

Statement

Without the written consent of Shenzhen ECGMAC Medical Electronics Co., Ltd (hereinafter referred to as ECGMAC), any materials contained in this manual shall not be photocopied, reproduced or translated into other languages. ECGMAC holds copyright in this manual and has the final interpretation right of this manual. ECGMAC owns the copyrights of this manual.

This manual contains proprietary materials protected by copyright law, including but not limited to technical secrets, patent information and other trade secrets. Users have the obligation of confidentiality and shall not disclose any content in this manual to any third party without written consent. The user's holding of this product manual does not mean that ECGMAC is willing to license the intellectual property rights contained in the product.

The pictures in this manual are for reference only, please be subject to the actual products. ECGMAC holds the rights to modify, update, and ultimately explain this manual.

It is reminded that the product shall be used strictly complying with this manual. User's operation failing to comply with this manual may result in malfunction or accident for which ECGMAC can not be held liable.

Responsibility of the Manufacturer

When using the equipment, please ensure that the electrical installation of the relevant room complies with national standards, and at the same time, follow the instructions in this manual. In the case that the following conditions are met, ECGMAC is responsible for any effect on safety, reliability and performance of the equipment.

- The equipment is used in accordance with the user manual;
- Assembly operations, extensions, re-adjustments, modifications or repairs are carried out by persons authorized by ECGMAC;
- The storage environment, working environment and electrical environment of the equipment meet the specifications;
- The label of the equipment is clearly identifiable, which can confirm that it was manufactured by ECGMAC;
- Damage not caused by non-human factors or force majeure factors such as typhoons and earthquakes.

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Shenzhen ECGMAC Medical Electronics Co., Ltd.

2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District,
518102 Shenzhen, China

Well Kang Ltd

Entreprise Hub, NW Business Complex, 1 Beraghmore Rd, Derry, BT548 8SE, Northern Ireland, UK

CONTENTS

1	Safety Guidance	52
1.1	Safety Warnings	52
1.2	Warnings and Cautions	53
1.3	Contraindication	53
1.4	Lithium Battery Care Warnings	54
1.5	List of Symbols	54
2	Product Introduction	55
2.1	Introduction	55
2.2	Composition	55
2.3	Product Details	56
2.4	Patient Cable Socket and Definition for Plug Pins	59
2.5	Patient Cable & Electrodes	60
3	Operation Preparations	60
3.1	Preparing the Equipment	60
3.2	Preparing the Patients	63
3.3	Connect the Patient Cable	63
3.4	Connect the Electrodes	64
4	Acquisition and Print ECG Report	65
4.1	Main Interface	65
4.2	Input Patient Information	66
4.3	Check the Quality of the ECG	69
4.4	ECG Acquisition	70
4.5	Reanalysis	71
4.6	Print ECG Report	71
5.	Setup	72
5.1	Print Setup	72
5.2	Analysis Setup	74
5.3	Communication Setup	75
5.4	Information Setup	76
5.5	Admin Configuration	76
5.6	Filter Setup	77
5.7	Time and Date Setup	79
5.8	Lead Selection Setup	79
5.9	Others	80
5.10	System Check	81
6	File Management	81
6.1	Upload Report	82
6.2	Data Exchange with External Storage	83
6.3	Check ECG Files on PC	84
6.4	Transmit Data to ECG Data Management Software	84
7.	Cleaning, Care and Maintenance	85
7.1	Overview	85
7.2	Cleaning	85
7.3	Disinfection	86
7.4	Care and Maintenance	86
8	Troubleshooting	88
8.1	ECG Waveform of Some Leads Cannot Be Displayed.	88
8.2	Keypad Not Work	88
8.3	AC Interference	88
8.4	EMG Interference	89
8.5	Baseline Drift	89
8.6	Save ECG Data without Printing	90
8.7	Paper-jam	90
9	Warranty and After-sales Service	90
9.1	Warranty	90
9.2	Manufacturer's Disclaimer	90
9.3	After-sale Service	91
9.4	Useful Life & Date of Manufacture	91
Appendix A	Package and Accessories	91
Appendix B	Technical Specifications	92
Appendix C	Key Components	95
Appendix D	List of Functional Accessories	95
Appendix E	EMC Information	95
Appendix F	Manufacturer Information	97
Appendix G	European Authorized Representative Information	97

1. Safety Guidance

1.1 Safety Warnings

- ⚠ The intended use is to acquire the ECG waveform from adults and children patients through electrodes.
- ⚠ The equipment should only be used in hospitals or medical institutions by doctors and trained healthcare professionals.
- ⚠ Please read this user manual carefully before operating.
- ⚠ The interpreted ECG with measurements and interpretive statements is offered to clinicians on an advisory basis only and require a doctor's signature to be effective.
- ⚠ The results given by the equipment should be examined based on the overall clinical condition of the patient, and they can not substitute for regular checking.
- ⚠ Please use the socket which has ground protection and ensure that the socket is properly grounded to avoid the risk of electric shock to the patient and operator.
- ⚠ Please ensure that the installation room has stable power supply system which is reliable and grounding.
- ⚠ When the power supply system is not complete and reliable, cut off the A/C power and directly use the internal D/C supply.
- ⚠ Do not use the equipment in the presence of flammable anesthetic mixtures with oxygen, hydrogen or other flammable agents.
- ⚠ Do not use in a medical high-pressure oxygen chamber, otherwise there is a danger of explosion.
- ⚠ Do not use this equipment in the presence of high static electricity or high voltage equipment which may generate sparks.
- ⚠ Only qualified service engineers can install this equipment, and only service engineers authorized by the manufacturer can open the shell.
- ⚠ Auxiliary equipments connected to the analog and digital interfaces must be certified according to the respective IEC/EN standards (e.g. IEC/EN 60950 for data processing equipment and IEC/EN 60601-1 for medical equipment). Furthermore all configurations shall comply with the valid version of the standard IEC/EN 60601-1-1. Therefore anybody, who connects additional equipment to the signal input or output connector to configure a medical system, must make sure that it complies with the requirements of the valid version of the system standard IEC/EN 60601-1-1. If in doubt, consult our technical service department or your local distributor.
- ⚠ When defibrillator is used simultaneously with the equipment, the operator should not touch patient, bed, table or equipment. All the electrodes (whether connected to the patient or not) and the patient do not need to be grounded. When the instrument is operated simultaneously with defibrillator or other electrical stimulation equipment, it is recommended to use disposable plate chest electrodes to avoid skin burns by metal electrodes.
- ⚠ If multiple instruments are connected to a patient, the sum of the leakage currents may exceed the limits given in the IEC/EN 60601-1 and may pose a safety hazard. Do not connect any equipment or accessories that are not approved by the manufacturer or that are not approved in IEC/EN 60601-1-1 to the instrument. The operation or use of non-approved equipment or accessories with instrument is not tested or supported. In this condition, the safety of instrument operation can not be guaranteed.
- ⚠ When the patient is equipped with pacemaker, it may influence the accuracy and results of the ECG examination. It is recommended that the doctor combines the waveform to

- ⚠ make the diagnosis. The existence of pacemaker will also increase the potential danger. In this case, pay special attention to safety when recording the ECG, proper measures must be taken to ensure that the leakage current is at a safe level.
- ⚠ To prevent burns, the contact points of the electrosurgical knife should be far from the electrode. The resistance between the electrosurgical knife and the patient's body should be as small as possible, and extra care should be taken. Plate electrodes can be used when necessary. Featured in their large contact area, high-frequency current density can be limited to an acceptable range.
- ⚠ Do not touch the patient, bed, table or the equipment while using the ECG together with a defibrillator or pacemaker.
- ⚠ Only the patient cable and other accessories supplied by the manufacturer can be used. Or else, the performance and electric shock protection can not be guaranteed.
- ⚠ Make sure that all electrodes are connected to the patient correctly before operation. Ensure that the conductive parts of electrodes and associated connectors, including neutral electrodes, do not touch the earth or any other conducting objects.
- ⚠ Don't touch the SIP/SOP and the patient simultaneously.

1.2 Warnings and Cautions

- ❗ The equipment should be placed on a flat surface avoiding excessive vibration and shock when moving.
- ❗ Avoid liquid splash and excessive temperature. The temperature must be kept between 5 °C and 40 °C during operation.
- ❗ Do not use the equipment in a dusty environment with bad ventilation or in the presence of corrosive substance.
- ❗ Make sure that there is no intense electromagnetic interference source around the equipment, such as radio transmitters or mobile phones etc. Attention: large medical electrical equipment such as electrosurgical equipment, radiological equipment and magnetic resonance imaging equipment etc. is likely to bring electromagnetic interference.
- ❗ Before using, the equipment, patient cables, and electrodes must be checked to seek any damage that may affect patient safety. If obvious damage or aging is found, replace the part.
- ❗ The frequency and voltage of the AC power supply shall be in accordance with the requirements and have sufficient current capacity.
- ❗ It is recommended to use the power adapter supplied by ECGMAC to avoid affecting the performance of the equipment and causing damages.
- ❗ The equipment should be placed in a quiet and comfortable environment.
- ❗ The device and accessories are to be disposed of according to local regulations after their useful lives. Alternatively, they can be returned to the dealer or the manufacturer for recycling or proper disposal.
- ❗ If there is any accident during the use, please stop to use it.
- ❗ To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Do not place ME equipment where it is difficult to operate the disconnect device.

1.3 Warnings and Cautions

None

NE 1.4 Lithium Battery Care Warnings

- ⚠ Improper operation may cause the internal li-ion battery (hereinafter called battery) to be hot, ignited or exploded, and it may lead to the decrease of the battery capacity. It is necessary to read the user manual carefully and pay more attention to warning messages.
- ⚠ DANGER OF EXPLOSION -- Do not reverse the anode and the cathode when installing the battery.
- ⚠ Do not use the battery near fire or in an environment where the temperature exceeds 60°C. Do not heat or splash the battery or throw it into fire or water.
- ⚠ Do not destroy the battery; do not pierce battery with a sharp object such as a needle; do not hit with a hammer, step on or throw or drop to cause strong shock; do not disassemble or modify the battery, otherwise the battery will heat, smoke, deform or burn, causing danger.
- ⚠ When leakage or foul smell is found, stop using the battery immediately. If your skin or cloth contacts with the leakage liquid, cleanse it with clean water at once. If the leakage liquid splashes into your eyes, do not wipe them. Irrigate them with clean water first and go to see a doctor immediately.
- ⚠ Stop using the battery when it reaches the end of its useful life, or if you notice a strange smell, deformation, discoloration, or distortion.
- ⚠ Properly dispose of or recycle the depleted battery according to local regulations.
- ⚠ Only when the device is off can the battery be installed or removed.
- ⚠ Remove the battery from the equipment when the equipment is not used for a long time.
- ⚠ If the battery is stored alone and not used for a long time, we recommend that the battery should be charged at least once every 6 months to prevent over discharge.
- ⚠ Keep the battery out of the reach of children.
- ⚠ Only authorized engineers can remove the battery box cover and replace the battery; and must use the same type of rechargeable lithium battery provided by our company.

1.5 List of Symbols

Symbols	Function	Symbols	Function
	USB socket		Alternating current
	DC Power Supply socket		Serial Number
	Manufacturer		Date and country of manufacture
	Recycle		Transportation and storage temperature range
	Avoid direct sunlight		Keep dry
	Tier limitation		DEFIBRILLATION-PROOF TYPE CF APPLIED PART
	Battery indicator		Note (general warning): the information you should know to avoid possible damage to patients or operators.
	Caution signs: Information you should know about how to avoid possible damage to your equipment.		Equipotentiality
	Humidity range of transportation and storage		Atmospheric pressure range for transportation and storage

	Fragile		Batch Number
	CE		Follow instructions for use
	European Authorized Representative		Recycling electronic equipment
	Medical device		Model number
	Importer		Distributor

2. Product Introduction

2.1 Introduction

Multi-channel ECG (EM-301&EM-601 series) acquire ECG signals from adults or children through the electrodes attached to human body and the ECG recorded by the equipment provide an important basis for analyzing and diagnosing heart disease. ECG signals of 12 leads are recorded simultaneously, and 12-channel waves are displayed and recorded simultaneously. The 12-channel ECG supports automatic measurement and automatic diagnosis functions. However the ECG with measurements and interpretation is offered to clinicians on an advisory basis only. Multi-channel ECG (EM-301&EM-601 series) includes EM-301, EM-301A, EM-301B, EM-601, EM-601A and EM-601B.

Scope of application: It is used by medical institutions to extract the human body's ECG wave for clinical diagnosis and research. The equipments intended to be used only in hospitals or healthcare facilities by doctors and trained healthcare professionals.

2.2 Composition


This equipment consists of the main unit, ECG data management software (optional), patient cable, grounding wire, power adapter, power cord, limb electrodes and chest electrodes. The ECG data management software is composed of ECG data replay, patient record management, report print and system setting modules. The software release version is V1, and it is loaded in the CD.

2.2.1 Product Models And Specifications

Model	Channel	Lead	Storage rate	Sampling rate	Display	Consumption (VA)	Dimension(mm)	Others
EM-301	3	12	1000sps/channel	2000sps/channel	7-inch color LCD display	≤40	310 244 65	12bit, with LAN Port
EM-301B	3	12	1000sps/channel	2000sps/channel	7-inch color LCD display		310 244 65	12bit, without LAN Port
EM-301A	3	12	1000sps/channel	32000sps/channel	7-inch color LCD display		310 244 65	12bit, with LAN Port
EM-601	6	12	1000sps/channel	2000sps/channel	7-inch color LCD display	≤40	310 244 65	12bit, with LAN Port
EM-601B	6	12	1000sps/channel	2000sps/channel	7-inch color LCD display		310 244 65	12bit, without LAN Port
EM-601A	6	12	1000sps/channel	32000sps/channel	7-inch color LCD display		310 244 65	12bit, with LAN Port

2.3 Product Details


2.3.1 Top Panel

Indicator	Name	Function
A	Recorder	Install recorder paper and print ECG waveform
B	LCD Screen	It displays the operation menu and supports 90-degree flipping
C	Keyboard	Selecting menu functions on the screen; Input numbers, characters and words
D	Paper casing switch	Press down the bottom can open the recorder case
E	Power Indicator	AC : Alternating current indicator Charge : Battery recharging indicator DC : If AC power is not connected, this indicator represents the DC power status. If the battery is being charged, the indicator is on alternately.

- ❗ Do not place heavy objects on the LCD screen or hit the LCD screen, otherwise, it will be damaged.
- ❗ When the equipment is not in use, put the display back to its original position to prevent accidental damage.


2.3.2 Bottom Panel

Indicators	Name	Function
A	Battery Compartment	Compartment for battery
B	Battery Label	The battery label indicates the rated voltage and the rated capacity of the rechargeable lithium battery.
C	Heat Emission Hole	Path for internal heat emission
D	Product Label	The product label indicates the product information.
The calibration output voltage and the calibration capacity of the rechargeable lithium battery pack are as follows : Calibration output voltage:14.4V; Calibration capacity : 2600mAh		


❗ If the lithium battery is not used for a long time (more than two to three months), the user should first charge the battery when using it again.

2.3.3 Rear Panel


Indicator	Name	Function
A	Power Supply Socket	Connected to power adapter
B	Heat Emission Hole	Path for internal heat emission
C	Handle	Portable

2.3.4 Right Panel

Indicator	Name	Function
A	Patient Cable Socket	Connect to the patient cable
C	SD Card Slot	Read the SD card
D	USB Socket	Standard USB socket, connecting to a U disk
E	USB Socket	Master USB socket, connecting to a PC to realize data exchange
F	Network cable port (optional)	Standard network cable port, connecting to a PC
G	Potential Equalization Conductor	When the equipotential grounding wire is needed to ensure the safety of electricity, the potential equalization conductor shall be connected with other grounding wires or installed grounding wires

2.3.5 Keyboard Function


Indicator	Name	Function
A	Alphanumeric keyboard	Input digit, letter and symbols
B	ESC	Cancel operation
C	Function Key	Select menu functions on the screen
D	Delete	Delete characters

E	Input method	Choose the input methods: English / Numbers
F	Arrow Keys	Moving the cursor (Up, Down, Left, Right)
G	COPY	Copy the last ECG signals when the system works in auto model
H	START/STOP	Start/Stop printing reports
I	Reset	A large polarization voltage may cause baseline drift. On the main screen, pressing the ESC key can decrease the polarization voltage and draw the baseline to zero quickly.
J	MODE	Press this key to select a working mode among the manual, auto,rhythm.
K	Lead Selection	In the manual mode, increase, decrease or switch patient cables
L	Power On/Off	Power-on/Power-off
M	SLEEP/WAKE UP	Rest/wake up the equipment

2.4 Patient Cable Socket and Definition for Plug Pins

2.4.1 Patient Cable Socket


: Applied part of type CF with defibrillator proof

2.4.2 Definitions of Corresponding Pins

Pin	Signal	Pin	Signal	Pin	Signal
1	C2(Input)	6	SH	11	F(Input)
2	C3(Input)	7	NC	12	NC
3	C4(Input)	8	NC	13	C1(Input)
4	C5(Input)	9	R(Input)	14	NC
5	C6(Input)	10	L(Input)	15	N or RF(Input)


2.5 Patient Cable & Electrodes

2.5.1 Patient Cable


The patient cable includes main cable and lead wires; which can be connected to electrodes according to the colors and identifiers. There are 6 chest lead wires and 4 limb lead wires.

2.5.2 Chest Electrode

2.5.3 Limb Electrode (Clamp type) :

3. Operation Preparations

3.1 Preparing the Equipment

Make sure that there is no intense electromagnetic interference source around the equipment, such as electrosurgical equipment, ultrasound instruments, radiographs, etc. Turn off the equipment if necessary.

The temperature must be kept between 5 °C and 40 °C during operation while humidity less than or equal to 80%.

3.1.1 Install Battery

Figure A

Figure B


Figure C


Figure D

Installation Method:

- 1) Turn the ECG machine upside down, rotate the fixing screw of batteries by anti-clockwise, and open the battery cover (as shown in Figure A above);
- 2) Take out the battery from the accessories box, insert the battery plug into the socket (the direction must be correct), and check if the wire array is consistent with Figure B;
- 3) Put the battery back into the battery compartment (as shown in Figure C above);
- 4) Close the battery cover and screw it clockwise (as shown in Figure D above).

3.1.2 Install Recorder Paper

Roll Paper : EM-601\EM-601A\EM-601B:110mm width
EM-301\EM-301A\EM-301B:80mm width

- 1) Press the button indicated in Figure A to open the recorder case.

Figure A




Figure B


- 2) Take out the paper roller, as shown in the figure below and then put the paper through the roller with the grid side of the paper facing the thermal print head;

3) As shown below, place the thermal paper and the roller gently in the recorder. Pull about 2 cm of paper out, and shut the recorder casing.

Loading Rolled Thermal Paper


Loading Rolled Thermal Paper

- ❗ When the recorder paper runs out or is not loaded, the hint message No Paper will appear on the screen. Then the user should load or replace the recorder paper immediately.
- ❗ Make sure that the recorder paper is installed in the center of the recorder, and the paper edge is parallel with the casing edge in the direction of advancing paper, in order to avoid paper deviation or damage to the paper edge.

3.1.3 Power Supply

Please check whether the power cord is connected to the unit well. The grounded three-slot outlet should be used. One end of the yellow equipotential grounding wire is plugged into the grounding hole of the ECG, and the other end is connected to the external ground wire.

⚠ Caution

Connect one end of the equipotential grounding wire to the grounding hole of the ECG and the other end to the ground to enhance the reliability of the grounding. Do not use water pipes or other pipes as the earth wire, otherwise the first-level safety measures of the equipment are invalid and patients may be at risk of electric shock.

3.1.3.1 AC Power On / Off

Connect the power cord, power adapter and equipotential grounding wire. When the

keyboard AC indicator light is on, press the «power»  key on the keyboard (about 3 seconds) and then release it, and then the instrument is turned on.

Press the «Power»  key (about 3 seconds) on the keyboard when the instrument is on. The content on the LCD screen disappears and enters shutdown state. Then, unplug the AC power cord and equipotential grounding wire.

3.1.3.2 DC Power On / Off

If the power cord is not connected, a built-in rechargeable lithium battery would be used, and DC light indicator is lit. Press the on the control panel directly to turn on/off the unit.

Battery charging :

When the instrument is connected to AC power supply, if the built-in battery is insufficient, the battery will be charged automatically, and the charging indicator light will be on when charging.

❗ Please follow the above instructions to turn on and off the machine, otherwise the instrument might not working properly.

❗ When the screen displays “The system is shutting down”, release the “Power” key.

3.2 Preparing The Patients

1. Before attaching the electrodes, please communicate with the patient and explain the procedure :

1) Before checking, patient should take a rest, keep clam and be relaxed :

2) Patient should not smoke, drink tea, coffee or alcohol before the check :

3) Patient should wear loosely if possible to facilitate the check :

4) Clean the skin where the electrode is placed;

5) Shave hair where the electrode is placed if necessary. Excessive hair prevents a good connection.

2. When possible, prepare the patient in a quiet room or area where others can't see the patient. If there is someone else in the room, please prepare a curtain when checking the patient.


3. Make sure that the patient is comfortable. The more relaxed the patient is, the less the ECG waveform will be affected.

4. Once the electrodes are attached to the patients, please inform the patient to keep normal breathing and keep still. Do not talk or chew.

3.3 Connect the Patient Cable

The Patient Cable will be used as applied parts specified.

As shown below, connect the patient cable to the patient cable socket on the right side of the main unit, and then screw it.

❗ The performance and electric shock protection can be guaranteed only if original patient cable and electrodes of the manufacturer are used.

❗ Patient cable socket can only be used to connect ECG signal, not for other purposes.


3.4 Connect the electrodes

It is important to connect the electrodes properly and ensure the electrodes are in good contact; new and old electrodes, reusable electrodes, and disposable electrodes cannot be mixed and used at the same time. The use of different types of electrodes will seriously affect the quality of the ECG recording. The electrode or cable plug must not touch other conductors, such as metal beds, etc., and the electrode should be replaced entirely when necessary.

3.4.1 Limb Electrodes Connection

Attach the limb electrodes to the electrode positions on body surface. First clean the areas where the electrodes will be placed with alcohol (75%) , then apply a small amount of gel to the cleaned skin, as shown in the figure :


R (RA) Right arm/Right deltoid;
L (LA) Left arm/Left deltoid;
RF (RL) Right leg/Upper leg as close to torso as possible;
F (LL) Left leg/Upper leg as close to torso as possible.

3.4.2 Chest Electrodes Connection

Use alcohol to clean the skin at the C1-C6 position of the chest, and then apply gel to the cleaned positions. Daub the round area of 25mm in diameter on each electrode site with gel evenly. Place a small amount of gel on the brim of chest electrode's metal cup. Place the electrode on the chest electrode site and squeeze the suction bulb. Unclench it and then the electrodes are affixed to the chest, the positions of chest electrodes on the body surface are as below:

C1 (V1) : Fourth intercostal space at the right border of the sternum;
C2 (V2) : Fourth intercostal space at the left border of the sternum;
C3 (V3) : Fifth rib between C2 and C4;
C4 (V4) : Fifth intercostal space on the left midclavicular line;
C5 (V5) : Left anterior axillary line at the horizontal level of C4;
C6 (V6) : Left midaxillary line at the horizontal level of C4.

❗ **Caution :** To avoid short circuit, gel should be placed separately and the chest electrodes should not be in contact with each other . If there is no gel, it can be replaced by alcohol (75%). Clean the electrode area on the chest surface with alcohol (75%). Connect the electrode to the limb, and make sure that the area where the electrodes placed is moderately moist. Do not replace the gel with saline to prevent corrosion of the electrodes.


3.4.3 Electrode Identification and Connection Color Table

Electrode Connectors	European		American Label	
	Identifier	Color Code	Identifier	Color Code
Right arm/Right deltoid	R	Red	RA	White
Left arm/Left deltoid	L	Yellow	LA	Black
Right leg/Upper leg as close to torso as possible	N or RF	Black	RL	Green
Left leg/Upper leg as close to torso as possible	F	Green	LL	Red
Chest 1	C1	Red	V1	Red
Chest 2	C2	Yellow	V2	Yellow
Chest 3	C3	Green	V3	Green
Chest 4	C4	Brown	V4	Blue
Chest 5	C5	Black	V5	Orange
Chest 6	C6	Purple	V6	Purple

4. Acquisition and Print ECG Report

4.1 Main Interface

Press “Power” for 1 second. After the device is turned on, the main interface pops up :

	Name	Explanation
A	Work Mode	By switching corresponding button, [Rhythm], [auto], [manual] mode can be chosen.
B	Baseline Filter	The value can be set in the Filter setup interface
C	EMG Filter	The value can be set in the Filter setup interface
D	AC Filter	The value can be set in the Filter setup interface
E	Time & Date	The value can be set in the Time & Date setup interface
F	Power supply indicator	Indicate the current power supply : AC () or DC(---)
G	Heart Rate	Real time Heart rate
H	Name	Patient name
I	ID	Patient ID
J	External storage	When connected with external storage, the SD or USB symbol will be displayed on the screen ;
K	Lead and wave-form display area	Display patient waveform
L	Gender	Patient Gender
M	Age	Patient Age
N	Paper Speed and Gain	The current paper speed and gain are displayed on the screen. When the gain set as 10/5 mm/mV, 20/10 mm/mV or Auto, the gain displayed is 10 mm/mV.
F1	Patient information	Press function key 【F1】 to open patient information interface
F2	Paper speed setup	Press function key 【F2】 , there are 6 kinds of paper speed to select, which are 5 mm/s, 6.25 mm/s, 10 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s
F3	Sensitivity Setup	Press function key 【F3】 , there are 9 sensitivities to select, which are 1.25 mm/mV, 2.5 mm/mV, 5 mm/mV, 10 mm/mV, 20 mm/mV, 40 mm/mV, 10/5 mm/mV, 20/10 mm/mV, AGC.
F4	Freeze	Press function key 【F4】 to freeze current ECG waveform. However, the freeze function would not available if ECG acquisition time less than 10s and the ECG waveform can be frozen up to 300seconds.
F5	More	Press function key 【F5】 to go to the next page

When the main screen is displayed, press the function key 【F5】 “More” to enter the second page of the main interface, as shown in the figure below.

	Name	Explanation
F1	Lead	Press function key 【F1】 to enter the Lead interface
F2	Display	Press 【F2】 to enter the display interface
F3	File	Press 【F3】 to enter the File management interface
F4	Setup	Press function key 【F4】 to enter the setup interface
F5	More	Press function key 【F5】 to go back to the main interface

4.2 Input Patient Information

When the main interface is displayed, press 【F1】 to enter patient information interface :

Patient

Operate : Select Switch

ID:

Inpatient ID:

First Name:

Outpatient ID:

Last Name:

Department:

Age:

Weight: kg

Date of birth: / / (Y/M/D)

Height: cm

Gender: Male

BP(S/D): / mmHg

Checked by:

Pacing: ☐ Off ☐ On

Confirmed by:

en

Confirm F1

Cancel F2

Clear F3

Copy F4

Worklist F5

4.2.1 Input Patient Information Manually

After input the patient information, press 【F1】 to go back to the main interface.

Explanation
Patient information can by typed : ID, Name, Birth date, Gender, Admission No, Outpatient No, Technician, Physician, Weight, Height and Blood pressure Three ways of ID input : 1. Accumulation : ID number are added in order from 1. 2. Time : ID is generated based on the acquisition time. 3. Input by users. ID can be input via keyboard or barcode scanning. Please refer to 4.2.2 for details. The default way of ID input is via time collection and it can be changed in system setting. Attention : The ID, Name, Inpatient No, Outpatient No, Technician and Physician should less than or equal to 63 ASCII characters.

4.2.2 Input Patient Information by Scanning Barcode

4.2.2.1 Barcode Setting

Select barcode 【F2】 to enter the barcode setting interface in the admin configuration interface, as shown in the following figure:

The screenshot shows the 'Barcode' configuration window. It has a title bar 'Barcode' and a right-side control 'Operate : Select Switch'. The main area contains a table with two columns: 'Start Address' and 'End Address'. The rows are labeled 'Item', 'ID:', 'First Name:', 'Last Name:', 'Sex:', 'Year Of Birth:', 'Month Of Birth:', 'Day Of Birth:', 'Male Code:', and 'Female Code:'. The values entered are: ID: 1, First Name: 9, Last Name: 1, Sex: 8, Year Of Birth: 8, Month Of Birth: 5, Day Of Birth: 8, Male Code: 1, and Female Code: 2. At the bottom, there are two buttons: 'Confirm F1' and 'Cancel F2'.

Item	Start Address	End Address
ID:	1	8
First Name:	9	13
Last Name:	1	8
Sex:	8	8
Year Of Birth:	8	11
Month Of Birth:	5	6
Day Of Birth:	8	9
Male Code:	1	
Female Code:	2	

4.2.2.2 Scan Barcode

- Scanner is connected to ECG machine through USB port
- Correctly set the bar code of the ECG
- ID generation mode is set by user input
- Bar code scanner should be aimed at the bar code. If you hear “Drip” prompt, the scan is successful.
- Press the patient information 【F1】 to enter the patient information interface to check the scan results

Example: According to the above setting, the result of scanning barcode 9787040195835 is as follows:

The screenshot shows the 'Patient' information window. It has a title bar 'Patient' and a right-side control 'Operate : Select Switch'. The main area contains various input fields for patient data. The values entered are: ID: 97870401, Admission No.: (empty), First Name: 95835, Outpatient No.: (empty), Last Name: 97870401, Department: (empty), Age: (empty), Weight: (empty) kg, Date of Birth: 1958 / 04 / 19 (Y/M/D), Height: (empty) cm, Sex: Male, BP (S/D): (empty) / (empty) mmHg, Check Dr: (empty), and Diagnosis Dr: (empty). At the bottom, there are five buttons: 'Confirm F1', 'Cancel F2', 'Clear F3', 'Copy F4', and 'Worklist F5'. An 'En' button is also visible in the bottom right corner.

Caution : if the barcode is set correctly, but only ID can be scanned, please try to scan the barcode of the «enter» key in the manual with the scanner, and then try to scan the barcode again.

4.2.3 Querying Patient Information through Worklist

When the patient information query protocol is set to Worklist, select “Worklist [F5]” to enter the Worklist interface as shown below:

Worklist		1/595		
ID	Name	Sex	Age	Time
15212960	Patient, Name	Others		20200117110124
01812022	Patient, Name	Others		20200116150801
86543027	Patient, Name	Others		20200116150801
92448565	Patient, Name	Others		20200116110121
81199639	Patient, Name	Others		20200116110120
09165901	Patient, Name	Others		20200116110120
96596515	Patient, Name	Others		20200108143052
79504955	Patient, Name	Others		20200108143051
97194232	Patient, Name	Others		20200108143051
11765367	Patient, Name	Others		20200108143050
51290419	Patient, Name	Others		20200108143050
14499502	Patient, Name	Others		20200103165516
67932598	Patient, Name	Others		20200103165512
06559509	Patient, Name	Others		20191228172559
93169160	Patient, Name	Others		20191228171109
23227367	Patient, Name	Others		20191228171059
73402212	Patient, Name	Others		20191228170927
51372583	Patient, Name	Others		20191228170915
07468991	Patient, Name	Others		20191227172234
90885012	Patient, Name	Others		20191227172234
Confirm F1 Cancel F2 Update F3				

1/87 : It means that the current first piece of patient information is selected, and a total of 87 pieces of information are queried. The blue background indicates the currently selected information.

Confirm F1 : Select the current information to enter the patient information

Cancel F2 : Cancel the query results and enter patient information manually

Update F3 : Enter the update range setting interface of Worklist, as shown below

Update range

Operate : Select Switch

Update mode: ☒ Auto ☐ Manual

Update date: ☒ Today

☐ Tomorrow

☐ Yesterday

☐ Date range

2020 / 03 / 10

—

2020 / 03 / 10

◀ ▶

Confirm F1

Cancel F2

Date mode	Auto : Automatically query the server to obtain the patient’s information list according to the set conditions; Manual : The patient information list is only obtained when the [Confirm] button is pressed in the update interface.
Update date	Today : Query today’s patient information Tomorrow : Query tomorrow’s patient information Yesterday : Query yesterday’s patient information Date range : Query the patient information of the set date range
Confirm F1	Query the server information according to the selected update time range, and return to check the query results when the server is set correctly.
Cancel F2	Cancel query and return. No results

- ! Caution :**
1. When use worklist function , please select the ID generation method as user input in the Admin configuration.
 2. The Worklist server information must be set correctly in the communication setup of the system setting for the successful Worklist update.

4.2.4 Querying Patient Information through HL7

When the patient information query protocol is set to HL7, F5 of the patient information interface is HL7.

Input the ID, press the HL7 button, query the HL7 server information, return to the query result and enter the patient information interface.


When use barcode scanner to put ID, the equipment will automatically query the HL7 server information and return to the query result and enter the patient information interface.

- ! Caution :**
1. For accurate ID input, please select the ID generation method as user input in the Admin configuration interface.
 2. The HL7 server information must be set correctly in the communication setup interface for the successful HL7 update.

4.3 Check the Quality of the ECG

4.3.1 Lead-off Alarm

The operator checks the connection status of the lead at any time. If the lead falls off, the display will prompt the electrode lead off alarm, and the corresponding lead waveform is a straight line, as shown in the figure below: C6 lead off.

Please carefully check whether the connection between the corresponding electrodes and human body, the corresponding electrodes and lead wire, and the main unit is reliable. When the connection is reliable, the prompt message will disappear.

- ! Note :**
1. If lead-off alarm appears during the acquisition, please re-sample.
 2. The data acquired in the case of lead off may lead to information loss, leading to missed diagnosis or misdiagnosis.

ZE 4.3.2 Eliminate Waveform Interference

If there is interference in the waveform, please check the filter setting to make it most suitable for the current environment.

If the interference can not be eliminated after setting the filter, check the patient's electrode connection and skin condition, and confirm that whether the ambient temperature is appropriate.

If none of the above steps can eliminate the interference, please check if there is any equipment with strong interference in the surrounding environment.

4.3.3 Pacemaker Detection Setting

The default setting for pacemaker detection sensitivity is off.

If patient equipped with pacemaker, please set the pacemaker detection to low in the others of Setup.

When the ECG acquisition for patient with pacemaker detection done, it is recommended to set pacemaker detection to off.

4.4 ECG Acquisition

When waveform is stable, click the «Start / Stop» button to start recording ECG data. When the acquisition starts, the main interface displays the word «Recording».

Three work modes are available: Manual mode, AUTO mode and Rhythm mode. Users can press **【Manual】**, **【Auto】** or **【RHY】** button to select.

AUTO mode	Press the start/stop to start recording, and the recording will be finished automatically. In simple mode 12 * 1 format, the acquisition time can be set. Switch the display and print format through the [Auto] key. File will be saved if acquisition done.
Manual Mode	Users can select any lead group to record ECG data. Start and stop ECG acquisition by pressing Start/Stop Key manually. Press Manual key to change the display and print format under manual mode.
Rhythm mode	Users can select rhythm lead group need to record ECG data. Start recording by pressing Start/Stop Key manually, and acquisition will be stopped and start to print automatically when recording time reaches 60seconds. Press RHY key to change the display and print format under rhythm mode.

4.4.1 AUTO Mode

Press the **【Auto】** button to enter the auto mode. The auto mode records 12-lead ECG data. In the Auto Mode, there are 6 display formats in the display setting interface: 3 4, 3 4+1R, 3 4+3R, 6 2, 6 2+1R, 12 1.

Press the START/STOP key to start printing. The default print time is 10 seconds. The print needs to be stopped or ended by pressing the **【START / STOP】** key.

4.4.2 Manual Mode

1) Press the **【Manual】** key to enter the manual work mode. By repeatedly pressing the [Manual] key, «Manual 1», «Manual 2», «Manual 3», «Manual 6», and «Manual 12» can be switched. Press **【◀】** and **【▶】** to switch the lead.

Manual1	Single channel mode. Record ECG data of single lead
Manual2	Two-channel mode. Record ECG data of two leads simultaneously
Manual3	Three-channel mode. Record ECG data of three leads simultaneously
Manual6	Six-channel mode. Record ECG data of six leads simultaneously

EM-601A, EM-601 and EM-601B support Manual1, Manual2, Manual3 and Manual6. EM-301A, EM-301 and EM-301B support Manual1, Manual2 and Manual3.

Press **【START / STOP】** to enable printing. The print duration is controlled by user.

Press **【START / STOP】** again to stop printing.

4.4.3 Rhythm Mode


Press the **【RHY】** to select the rhythm work mode. The equipment acquires data for the specified rhythm lead. The acquisition time is 60 seconds, which has single lead rhythm and 3-lead rhythms. Press **【◀】** and **【▶】** to switch the lead.

Rhythm 1	Record the ECG data of selected single lead.
Rhythm 3	Records the ECG data of selected three leads simultaneously

After the waveform is stable, press **【START / STOP】** to start printing. The print duration is 60s. Press **【START / STOP】** again to stop print.

4.5 Reanalysis

Press F4 to enter «Setup» → «Analysis Setup», and enable «Reanalysis» function.


Confirm F1	Save the parameter modification and exit, and print the average template and modified parameters and interpretations at the same time.
Cancel F2	Abandon recording and printing.
Diagnosis F3	Display the interpretation of waveform.
Wave F4	Enter replay interface.
Enter	Switch the selected calibration line.
*	Set the move step of the calibration line. The available move steps are 2ms, 4ms, 8ms and 16ms.
<>	Move the calibration line left and right according to the set steps

4.6 Print ECG Report

This equipment supports two print methods: thermal printer and USB printer. You can select the appropriate printing method in the print setup interface.

thermal printer	Install thermal print paper and print the ECG report on the thermal paper
USB printer	Connect the printer and ECG machine through a USB cable, and print the ECG report on A4 paper.

Under Manual Mode, Rhythm mode, only thermal printer is supported, while under Auto Mode, both thermal printer and USB printer are supported.


3 4+3R

3 4+1R

5. Setup

Press **【F4】** to enter the Setup interface, as shown below:


In the Setup interface, press the corresponding number key to enter the corresponding setting interface. For example, press **【1】** to enter the print setup interface.

Restore factory settings/Default F2: Restore the related parameters of the machine to the factory default state.


5.1 Print Setup

Press button **【1】** to enter print setup interface in the system setting interface. Use **【▲】** and **【▼】** button to select Printer. The interface is shown as the figure below:

5.1.1 Thermal Printer


Press button **【1】** to enter thermal printer setup interface in the print setup interface. The interface is shown as the figure below:

ECG Grids Print	On, Off
Sample mode	Pre-sample: 10s ECG data sampled before pressing the START/STOP key will be printed out. Real-time Sample: ECG data sampled after pressing the START/STOP key will be printed out. The data will be recorded as soon as the user pressing the "START/STOP" button. Trigger Mode: After pressing the START/STOP key, if arrhythmia ECG data is detected during the acquisition process, the print will be triggered automatically. Otherwise, nothing will come out. Periodic Sample: Sample mode is in periodic for the same data. For example, if Interval is set to 2 min, and Duration is set to 24 min, after pressing the START/STOP key, the print will be performed every two minutes and come to 12 times.
Duration	Set the duration of periodic sample
Interval	Set the interval of periodic sample

5.1.2 USB Printer


Press button **【2】** in the print setup interface to enter the USB printer setup interface, as shown below:

ECG Grids Print	ECG grid print Off: Do not print the ECG grid Dots: The ECG grid is a grid of dots Line: The ECG grid is a net of lines
Print Direction	Under USB Printer Mode, print report direction can be horizontal or vertical.
USB printer model	When the USB cable and printer are connected, it will automatically check the model of the USB printer. If the ECG does not support the current printer model, a hint will pop out, and the report cannot be printed. 1. The USB printer will not immediately start the print at the end of the acquisition. Please be patient. 2. If the printer still cannot start the print after 30s, please check whether the connection between the ECG and the printer is valid, and whether the ECG supports this type of printer.

5.2 Analysis Setup

In the System Setting interface, press button **【2】** to enter the Analysis Setup, as shown below:


ECG Measurement	Measurement information will be showed in the printout or files when the ECG measurement is enabled. The measurement information includes: HR---Heart Rate PR interval QRS complex duration QT/QTc P/QRS/T RV5/SV1 RV5+SV1
-----------------	---

ECG Classification	ECG classification will be showed in the printout files when the ECG measurement enabled such as "Normal ECG".
ECG Diagnosis	Interpretation will be showed in the printout ECG files when ECG diagnosis is enabled.
Print detailed measurement parameters	On: print or save the current ECG detailed parameters of template of each lead displayed in picture. Off: the detailed parameters of template are not displayed.
Reanalysis	On: Reanalysis interface will be showed in auto mode once ECG acquisition is done if Reanalysis is enabled. Off: Acquire 10S data in auto mode and the equipment will directly print ECG analysis interpretations.

5.3 Communication Setup

In the Setup interface, press button 【3】 to enter the Communication Setup, as shown below:

The corresponding setting interface can be accessed through the digital key. The premise of successful communication setting is that the ECG must be connected to the Internet.

IP setup	In the IP setup interface, the default one is Auto obtained IP. If the user needs a fixed IP, please select a static IP and input the corresponding value. Note: When setting a static IP address, the user must ensure that the set static IP is not consistent with the IP of other devices in the LAN.
FTP setup	To realize the FTP upload function, the FTP setup of ECG machine must be set according to the FTP server information.
Shared directory settings	To achieve the shared directory upload function, the shared directory setup in the ECG must be set. The corresponding information of shared directory setup must be input according to shared directory of the computer.
DICOM setup	To realize the DICOM upload function, the DICOM setup in the ECG must be set. DICOM setup input the corresponding information based on the DICOM server.
HL7 setup	To access the HL7 server, the HL7 server of the ECG must be set. The specific settings are based on the HL7 server.
WORKLIST setup	To access the WORKLIST server, the WORKLIST server of the ECG must be set, and the specific settings are based on the WORKLIST server.

5.4 Information Setup

In the Setup interface, press button **【4】** to enter the Information setup, as shown below:


The screenshot shows the 'Information setup' screen. At the top, it says 'Information setup' and 'Operate : Select Switch'. There are two input fields: 'Hospital name:' and 'Device name:'. Below these fields are two buttons: 'Confirm F1' and 'Cancel F2'. There is also a small 'en' button in the bottom right corner.

In the Information setup, the hospital name and device name can be set. The set hospital name will be displayed on the printed report.

5.5 Admin Configuration

In the Setup interface, press button **【5】** on the keyboard and then the password input prompt box will pop out. Input the password (default password 123456) and then press **【F1】** to enter the admin configurations, as shown below:

The screenshot shows the 'Admin configuration' screen. At the top, it says 'Admin configuration' and 'Operate : Select Switch'. The screen contains various settings with radio buttons:

- ID generated by: ☐ Acc. ☒ Time ☐ User input
- Upload: ☐ Auto ☒ Manual
- Upload Method: ☐ FTP ☒ Shared Directory ☐ DICOM
- File Format: ☐ FDA-XML ☒ EM-XML ☐ SCP ☐ BKG ☐ DICOM
- Picture Type: ☐ Don't save ☒ PDF ☐ JPG ☐ BMP
- Picture Direction: ☐ Horizontal ☒ Vertical
- Picture ECG Grid: ☐ Off ☐ Dot ☒ Line
- Picture Grid Color: ☐ Gray ☐ Blue ☒ Red
- Patient Query Protocol: ☐ Off ☐ HL7 ☒ Worklist
- Data Mode: ☐ ECG mode ☒ Demo ☐ Calibration

At the bottom, there are five buttons: 'Language F1', 'Logo F2', 'Update F3', 'Barcode F4', and 'More F5'.

Item	Sub-project name	Description
ID generated by	Acc.	ID accumulated from 1 and not editable
	Time	ID generated automatically according to the check time. The ID length is 12 digits and it is not editable
	User input	ID input by user
Upload	Auto	File will be upload automatically
	Manual	The user selects the file to be uploaded and clicks the upload button to complete the upload
Upload protocol	FTP	To upload through FTP protocol, ECG and FTP server must be in the same LAN
	Shard Directory	To upload through the shared directory protocol, the ECG and the computer of shared directory must be in the same LAN
	DICOM	To upload through DICOM protocol
File format	FDA-XML/ EM-XML/ SCP/ BKG/ DICOM	
picture type	Don't save/ PDF/JPG/ BMP	
Picture direction	horizontal	Generate horizontal reports
	Vertical	Generate vertical reports
Picture ECG Grid	off	Generate report without ECG grid
	Dot	The ECG report generated has ECG grid formed by dot
	Line	The ECG report generated has ECG grid formed by lines
Picture Grid color	gray	The color of grid is gray / blue / red
	blue	
	red	
Patient query protocol	WORKLIST	Access to the Worklist server to check the patient protocol
	HL7	Access to the HL7 server to check the patient protocol
Data mode	ECG mode	Real time ECG waveform is displayed in the interface
	Demo	DEMO ECG waveform is displayed in the interface (heart rate of 80 bpm), and the word "demo" is displayed on the screen. In this mode, real-time signals cannot be obtained.
	Calibration	Calibration waveform is displayed in the interface (1mv square wave). In this mode, real-time signals can not be obtained.

Press button **【F5】** "More" to enter the second page of Admin Configuration. The second page and functions are as follows:

Confirm	Save the settings and exit the Admin Configuration
Cancel	Cancel the settings and exit the Admin Configuration
Change P	Change Admin Configuration password
Sys check	Enter the system diagnosis to self-check the hardware of the ECG, see 5.10 for details
More F5	Function button for switching to the next page of the Admin Configuration

Admin configuration Operate :

ID generated by: ☐ Acc. ☒ Time ☐ User input

Upload: ☐ Auto ☐ Manual

Upload Method: ☐ FTP ☐ Shared Directory ☒ DICOM

File Format: ☐ FDA-XML ☐ EM-XML ☐ SCP ☐ BKG ☐ DICOM

Picture Type: ☐ Don't save ☐ PDF ☐ JPG ☐ BMP

Picture Direction: ☐ Horizontal ☐ Vertical

Picture ECG Grid: ☐ Off ☐ Dot ☐ Line

Picture Grid Color: ☐ Gray ☐ Blue ☐ Red

Patient Query Protocol: ☐ Off ☐ HL7 ☐ Worklist

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Language F1	Enter the language switching interface, the user can set the machine language.
Logo F2	Press 【F2】 and input the correct password to enter the LOGO setup. The set LOGO will be displayed on the printed report. The maximum input length of LOGO setup is 16 characters.
Upgrade F3	The machine will automatically turn off when upgrading. Note: upgrading needs to be carried out under the guidance of professionals!
Barcode F4	In the barcode setting interface, the users can set the start address and end address of each parameter of the patient information. External barcode scanner needed to be connected and the barcode needed to be set.
More F5	Function button for switching to the next page of the Admin Configuration.

5.6 Filter Setup

In the Setup interface, press button **【6】** to enter the filter setup. At the same time, the user can also press the filter shortcut key on the main interface to enter the filter setup.

Filter Operate :

AC Filter: ☐ Off ☒ 50Hz ☐ 60Hz

Baseline Filter: ☐ Off ☐ 0.05Hz ☐ 0.16Hz ☐ 0.25Hz

☐ 0.32Hz ☐ 0.5Hz ☐ 0.67Hz

EMG Filter: ☐ Off ☐ 20Hz ☐ 25Hz ☒ 30Hz

☐ 35Hz ☐ 40Hz ☐ 45Hz

Lowpass Filter: ☐ Off ☒ 75Hz ☐ 100Hz ☐ 150Hz

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Item	Description
AC Filter	OFF, 50Hz, 60Hz
Baseline Filter	OFF, 0.05 Hz, 0.16 Hz, 0.25 Hz, 0.32 Hz, 0.5 Hz, 0.67 Hz
EMG Filter	OFF, 20 Hz, 25 Hz, 30 Hz, 35 Hz, 40 Hz, 45 Hz
Low pass Filter	OFF, 75 Hz, 100 Hz, 150 Hz

5.7 Time and Date Setup

In the Setup interface, press button 【7】 to enter time and date setup interface, as shown below:

Date format	yyyy-mm-dd/ mm-dd-yyyy / dd-mm-yyyy;
Time format	24-hour 12-hour(there are AM and PM for 12-hour format)
Date	Set the current specific date of the machine
Time	Set the current specific time of the machine

5.8 Lead Selection Setup

In the Setup interface, press button 【8】 to enter the Lead Selection Setup, the users can set the rhythm lead and lead system, as shown below:

There are 8 lead setup systems available for selection. Users can select one according to their needs.

Name	Lead system
Standard lead system	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
Posterior wall of the lead system	I, II, III, aVR, aVL, aVF, V1, V2, V3, V7, V8, V9
R-pectoral Lead system	I, II, III, aVR, aVL, aVF, V1, V2, V3R, V4R, V5R, V6R
Right chest lead system	I, II, III, aVR, aVL, aVF, V3R, V4R, V5R, V7, V8, V9
Intercostal lead system	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
Next intercostal lead system	I, II, III, aVR, aVL, aVF, V.1, V.2, V.3, V.4, V.5, V.6
Cabrera Lead System	aVL, I, aVR, II, aVF, III, V1, V2, V3, V4, V5, V6
User defined Lead	User defined

Single lead rhythm: Select one lead of the 12 leads as the rhythm lead. In the Auto mode, 1R is chosen in the format of 3 4+1R or 6 2+1R.

Three-lead rhythm: Choose 3 leads of 12 leads as the rhythm leads. In the Auto mode, 3R are chosen in the format of 3 4+3R.

5.9 Others

In the Setup interface, Press button **【9】** to enter others setup, as shown below:

Others Operate : < Select > Switch

Paced Detection Sensitivity: ☐ Off ☒ Low ☐ High

ECG color (display): ☐ Green ☐ Yellow ☐ White

Paper speed display: ☐ Common ☐ All

Gain display: ☐ Common ☐ All

Automatic Gain: ☐ Off ☐ On

Auto Update U-disk File Mapping: ☐ Off ☐ On

Keypad tone: ☐ Off ☐ On

Beep Volume: ☐ Off ☐ On

QRS Beep: ☐ Off ☐ On

Power-save: 0 min

LCD backlight off time: 0 min

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
Paced Detection Sensitivity	Off Low High
ECG color on Display	Green, Yellow, White
Paper speed display	Switch the paper speed in the main interface Normal : 12.5 mm/s, 25 mm/s, 50 mm/s All : 5 mm/s, 6.25 mm/s, 10 mm/s, 12.5 mm/s, 25 mm/s and 50 mm/s
Gain	Switch the gain in the main interface Normal : 5 mm/mV, 10 mm/mV, 20 mm/mV; All : 1.25 mm/mV, 2.5 mm/mV, 5 mm/mV, 10 mm/mV, 20 mm/mV, 40 mm/mV, 10/5 mm/mV, 20/10 mm/mV, AGC

Automatic gain	It is only valid for gain of 10 mm/mV. This function can be triggered only when the amplitude of input signal reach the automatic gain standard.
Auto update U-disk File Mapping	Turn on: When the master USB cable is connected to the ECG and the computer, the display file in the U-disk mapped by the computer-side ECG will be automatically updated after the acquisition is completed. Turn off: The display file in the mobile hard disk mapped by the ECG machine on the computer side remains unchanged. If it needs to be updated, please plug in the USB cable again. Note: this function will reduce efficiency of the system. Please use with caution.
Keypad tone	Whether the keyboard tone is turned on tone
Beep volume	The software prompts whether to turn on the beep volume
QRS Beep	QRS Beep can be set on and off here.
Power-save	If the machine is set aside without being operated, it will shut down within the set time and only numbers can be input here. For example, input 5, machine will automatically shut down if the users do not operate the machine for five minutes.
LCD backlight off time	If the machine is set aside without being operated, the display will be turned off within the set time and only numbers can be input here. For example, input 5, machine will automatically turn off the display if the users do not operate the machine for five minutes.

5.10 System Check

In the Admin configuration interface, select the sys check **【F4】** to enter the System check interface, users can test the machine in terms of keyboard, file system, battery, amplify board, net and printer here.

Press **【Enter】** to enter the test item. After the test, press the **【ESC】** to exit the current test item interface and return to the System check interface.


6. File Management

On the second page of the main interface, press File **【F3】** to enter the file management interface. Patient records can edited, printed, queried, previewed, deleted, and uploaded here.

File		5/1000						
ID	Name	Sex	Age	Check Time	File Format	Upload	Exp/Imp	
20200628092820	666666 666666	Female	555	20200628092820	DICOM			
1		Male		20200505114244	EM-XML			
1		Male		20200505114040	EM-XML			
20200323174952		Others		20200324121920	DICOM			
20200323174952		Others		20200324121818	DICOM			
20200323174952		Others		20200324121714	DICOM			
20200323174952		Others		20200324121611	DICOM			
20200323174952		Others		20200324121508	DICOM			
20200323174952		Others		20200324121406	DICOM			
20200323174952		Others		20200324121303	DICOM			
20200323174952		Others		20200324121200	DICOM			
20200323174952		Others		20200324121057	DICOM			
20200323174952		Others		20200324120955	DICOM			
20200323174952		Others		20200324120851	DICOM			
20200323174952		Others		20200324120749	DICOM			
20200323174952		Others		20200324120646	DICOM			
20200323174952		Others		20200324120543	DICOM			
20200323174952		Others		20200324120441	DICOM			
20200323174952		Others		20200324120338	DICOM			
20200323174952		Others		20200324120235	DICOM			
Transmit F1 Delete F2 Preview F3 Back F4 More F5								

F1Transmit	Files not uploaded, All files, Selected files  upload successful;  upload not successful
Delete F2	Delete "uploaded files", "all files", "selected files"
Preview F3	Press F3 to enter Preview
Back F4	Press F4 to exit File management
More F5	Press F5

File management interface : Press **【F5】** "More" to enter the second page. The second page buttons and functions are as follows:

Edit F1	Press F1 to enter the patient information interface. Users can modify the patient information by using the numeric keys and letter keys
Query F2	Press F2 to enter the query interface. In the query interface, users can search for files according to information such as ID, name, gender, age, file format, collection time, etc. When the time is on, only the file information within the set time is queried.
Export F3	Export files from ECG machine to external storage. SD card and U disk are supported Status of exported files: selected files / not exported / all files  It indicates that the file has been exported to external storage When exporting data, the users can choose to "retain" or "delete" the exported internal storage file
Import F4	Realize the function of importing files from external storage Select the storage where the imported file is located. SD card and U disk are supported  It indicates that the file is internally stored  It indicates that the file was imported from external storage Note: The import file format must be the format supported by the equipment, otherwise the import will fail.
More F5	Press F5 to show more function keys.

6.1 Upload Report

The report can be uploaded through FTP, shared directory, DICOM and other protocols.

6.1.1 Upload Setup

Press F4 in the main interface to enter «System Setting» → «Admin Configuration», and select the upload method and upload mode. Set the uploaded file format and image format at the same time.

When save the picture is set, the file and the corresponding picture report will be uploaded at the same time; while not save the picture is set, only the file will be uploaded.

! Note :

1. The upload method must correspond to the set file format, otherwise it cannot be uploaded.
2. To achieve DICOM upload, the file format must be set to DICOM, otherwise it cannot be uploaded.

6.1.2 Communication Setup

1. Press F4 in the main interface to enter “System Setting” → “Communication Setup” → “IP Address Setup” to connect to the network
2. Set the corresponding server information according to the upload method.

FTP: Enter «FTP Setting» to set FTP server information.

DICOM: Enter «DICOM Setting» to set DICOM server information

Shared Directory: Enter “Shared Directory Setup” to set shared directory information.

6.1.3 Upload

Upload is set to automatic: After the acquisition, ECG files and picture reports are automatically uploaded

Upload is set to manual: After the acquisition, go to «File Management», click «Upload», select the category of the uploaded file, and click «OK» to realize the upload.

! Note :

1. When there is upload error prompt, please check the equipment and network setting of server and related information.
2. In order to upload the file successfully, the file format and server settings should be correct. In addition, the status of server and the successful internet connection of the equipment should be ensured.

6.2 Data Exchange with External Storage

6.1.3 Upload

1. Insert U disk or SD card into ECG machine
2. Enter to «File Management» and click «Input»
3. Select the appropriate import media and click «OK»

! Note :

1. Only import file formats supported by the equipment.
2. Only ECG data file instead of picture can be imported.
3. When importing files, users must select the same properties as the inserted external storage. For example, select SD card to insert SD card, and you select U disk to insert U disk.
4. The system format of U disk or SD card must be FAT32, otherwise the files will not be imported.
5. When importing data, users must ensure that the external storage is effectively connected to the ECG machine, otherwise the import failure may occur.
6. Users must insert a USB flash disk or SD card before entering file management, otherwise the import may fail.

6.2.2 Data Export

1. Insert U disk or SD card into ECG machine.
2. Enter «File Management» interface and click «Export».
3. Select the appropriate export file type and click «OK».

! Attention :

1. The system format of U disk or SD card must be FAT32, otherwise the files cannot be exported.
2. When exporting files, users must choose the same attributes as the external storage. For example, select the SD card when insert the SD card, and select the U disk when insert the U disk.
3. Insert the U disk or SD card before entering the file management, otherwise the import may fail.

6.3 Check ECG Files on PC

When ECG machine connected to PC via USB cable, users can check the ECG files on PC.

! Note :

1. Disk mapping can only view file instead of editing.
2. Automatic U disk file mapping will reduce the efficiency of the ECG, please use it with caution.

6.4 Transmit Data to ECG Data Management Software

6.4.1 Transmit Files to ECG Data Management Software through FTP

1. Install the ECG data management software to the computer as an administrator.
2. Run the ECG data management software as an administrator.
3. ECG machine and PC are connected to the same network or connect the ECG network port with the computer network port through a network cable.
4. Press Setup-->Admin Configuration-->Communication setup to choose FTP upload as the upload method.
5. Input the FTP client information of ECG data management software in system setting-communication setup-FTP setup of ECG machine.
6. When the upload mode is selected to be auto, the file will be automatically uploaded after the acquisition; when the upload mode is selected to be manual, the users need to go to the file management interface to select the file to be uploaded.
7. After successfully upload, click the «Query» button on the database interface of the ECG data management software on the computer to check the uploaded file.

! Note :

1. In order to upload the file successfully, please make sure that the ECG data management software is working normally.
2. The ECG data management software is optional.
3. When the ECG machine and the computer of ECG data management software are not in the same network segment, please contact the manufacturer for technical support.
4. If the two ends of the network cable are used to connect the ECG network port and the computer network port, the first three segments of IP address are consistent with the computer, and the fourth segment is not the same. The default gateway and subnet mask are the same with the computer.

6.4.2 Import Files to ECG Data Management Software

After the ECG data of the ECG machine is imported into the computer installed with the ECG data management software, run the ECG data management software as an administrator. Click the **【Import】** button in the database interface of the ECG data management software to set the import path as the path of the ECG data. Select the ECG data to be imported and click confirm, then the ECG data is successfully imported to the ECG data management software.

7.Cleaning, Care and Maintenance

7.1 Overview

Please keep the equipment and its accessories free of dust. To prevent damage to the device, please follow these steps :

- ❗ Turn off the power before cleaning and disinfection. The mains supply must be switched off if it is in use and unplug the power cord and patient cables.
- ❗ Prevent the detergent from seeping into the main unit while cleaning. Do not immerse the unit or the patient cable into liquid under any circumstances.
- ❗ Do not clean the unit and accessories with abrasive fabric and avoid scratching the electrodes.
- ❗ After cleaning, remove detergent remainder with a clean dry cloth.
- ❗ This equipment and its accessories must be regularly maintained and checked (at least once a year).
- ❗ ECG machines are classified into Measuring Instrument, so users should send it to Official Measurement Administrative to test and certify according to national metrological calibration regulation of electrocardiograph and electroencephalogram machine every year.
- ❗ The instrument signal input/output connector (when needed) must be connected with Class I equipment which is EN 60601-1:2006 compliant, and the total leakage current should be tested to be available by users themselves.
- ❗ Only qualified service engineers can service and repair this equipment. When the equipment has malfunction, it should be clearly marked to avoid running with failure.
- ❗ The equipment can not be modified or altered in any way.

7.2 Cleaning

7.2.1 Cleaning the Main Unit

Use a soft cloth dipped in water, neutral detergent or alcohol (75%) to clean the external surface of the main unit after dry it. Wipe the socket, paper compartment and operation panel of ECG machine with dry cloth. Clean the display screen with a dry and soft cloth or a dry cloth dipped in neutral detergent. Avoid water entering the machine, which may cause malfunction.

7.2.2 Cleaning the Patient Cable and Electrodes

Clean the patient cable and electrodes with a soft cloth dipped in water, neutral detergent or alcohol (75%), and dry them with a dry cloth to ensure that the patient cable and electrodes are completely dry after cleaning.

- ❗ Note: if there is gel on the surface of patient cable or electrodes or electrode pins or metal parts are wet, inaccurate recording of ECG waveform may occur.

7.2.3 Cleaning the Print Head

Dirty and soiled thermal print head will deteriorate the printing definition. So it should be cleaned at least once a month regularly.

Open the recorder casing and remove the recorder paper. Wipe the print head gently with a clean soft cloth dampened in alcohol (75%). For stubborn stain, soak it with a little alcohol (75%) first and wipe it off with a clean soft cloth. After air dried, install the recorder paper and close the recorder casing.

7.2.4 Cleaning the Silicone Rubber Shaft in the Recorder Casing

The silicone rubber shaft should be flat, smooth, and smudge-free, otherwise the printing is not clear. Wipe in the longitudinal direction with a clean soft cloth damped in alcohol (75%). Close the recorder casing after the alcohol evaporate.

7.3 Disinfection

Before disinfection, clean the equipment first. Then wipe the surfaces of the unit and the patient cable with hospital standard disinfectant.

- ❗ Do not use high temperature, high pressure steam, or ionizing radiation for disinfection.
- ❗ Do not use chloric disinfectant such as chloride, sodium hypochlorite etc.

Turn off the power and wipe the surface of main unit, patient cable and electrodes with a clean and soft cloth damped in disinfectant (alcohol 75%). Dry the equipment with a clean dry cloth or air dry naturally.

7.4 Care and Maintenance

- ❗ Only qualified service engineers can repair this equipment. When the equipment has malfunction, it should be clearly marked to avoid running with failure.

7.4.1 Recharge and Replacement of Battery

7.4.1.1 Recharge of Battery

The ECG machine has a built-in rechargeable lithium battery and charge control circuit. When the battery is used for the first time, because of the loss of power in the process of storage and transportation, the battery is generally not sufficient. Before use, the battery should be charged in advance. As long as the AC power supply is connected, the battery can be charged. At this time, the AC power indicator and the battery charging indicator light are on at the same time, indicating that the battery is charging. When the battery is fully charged, the battery charging indicator is off. The charging time from the exhausted state to 90% of the battery capacity is about 7.5 hours. Under the condition of battery power supply, the continuous working time of the ECG machine should not be less than 1 hour.

7.4.1.2 Capacity Identification

The current capacity of the battery can be identified according to the battery symbol in the top right corner of the LCD screen, as follow:

-  Full capacity
-  Capacity is adequate.
-  Capacity is low and needed to be charged.
-  Capacity is insufficient and needed to be charged.
-  Capacity is exhausted and needed to be charged immediately.

7.4.1.3 Battery Replacement

When the useful life of the battery is over, or foul smell and leakage are found, please contact the manufacturer or the local distributor for replacement.

7.4.2 Recording Paper

It is recommended to use recording paper provided by the manufacturer. The usage of other paper may shorten the life of the thermal print head. The deteriorated print head may lead to illegible ECG reports and block the advance of the paper.

Tips for choose and Store Recording paper	1.Never use wax-coated, gray-black recorder paper to prevent the wax from sticking to the print head, which may cause malfunction or damage.
	2.Recording paper should be stored in a dry, dark and cool area, avoiding excessive temperature, humidity and sunshine. Do not place the recording paper under the fluorescent lamp for a long time to avoid affecting the recording effect.
	3.Make sure that there is no polyvinyl chloride or other chemicals in the storage environment, which will lead to color change of the paper.
	4.Do not overlap the recorder paper for a long time, or else the ECG reports may trans-print each other.
	5.Special attention should be paid to the size of the recorder paper. Recorder paper that does not meet the requirements may damage the thermal print head or the silicone rubber shaft.

7.4.3 Maintenance of the Main Unit and the Patient Cable

Besides the maintenance requirements recommended in this manual, please comply with local regulations of maintenance and measurement.

The following safety checks should be performed at least every 12 months by a qualified person who has adequate training, knowledge, and practical experience to perform these tests:

- Inspect the equipment and accessories for mechanical and functional damage:
- Inspect the safety related labels for legibility:
- Inspect the fuse to verify compliance with the rated current and circuit-breaking characteristics:
- Verify that the device functions properly as described in the instructions for use:
- According to EN 60601-1:2006 run the following test:

Protection earth resistance, Limit: 0.1 Ω

Earth leakage current, Limit: NC 500 μ A, SFC 1000 μ A

Enclosure leakage current, Limit: NC 100 μ A, SFC 500 μ A

patient leakage current, Limit: 10 μ A (CF Type)

patient auxiliary current, Limit: NC a. c. 10 μ A, d. c. 10 μ A; SFC a. c. 50 μ A, d. c. 50 μ A

patient leakage current under single fault condition with mains voltage on the applied part, Limit: 50 μ A (CF Type)

Only qualified person who has adequate training, knowledge can perform these tests and makes record; If the device is not functioning properly or fails any of the above tests, the equipment has to be repaired.

7.4.3.1 Main Unit

- Disconnect the power supply to power off the equipment;
- Wipe the main unit and its accessories. Put the dustproof coat on the main unit after use;
- Prevent any liquid from seeping into the equipment, otherwise the safety and performance of the equipment can not be guaranteed;
- Ask the maintenance department to check the equipment regularly.

7.4.3.2 Patient Cable

- Check whether the patient cables are in good contact according to the following table. The resistance between each lead from the electrode plug to the corresponding pin in the lead wire plug is less than 10 Ω .
- Note: The resistance of the patient cable with defibrillation-proof is approx.10 K Ω .

connector symbol	R	L	F	RF	C1	C2	C3	C4	C5	C6
connector pin number	9	10	11	14	12	1	2	3	4	5

- Integrity of the patient cable, including the main cable and lead wires, should be checked regularly. Make sure that it is conductible;
- Arrange the patient cable to avoid twisting, knotting or crooking in a closed angle while using it;
- Do not drag or twist the patient cable with excessive strength while using it. Hold the connector plug instead of the cable when connecting or disconnecting the patient cable;
- Store the lead wires in a big wheel or hang it to avoid twisting;
- Once damage or aging of the patient cable is found, replace it with a new one immediately.

7.4.3.3 Electrodes

- Electrodes must be cleansed after use and make sure there is no remainder gel on them;
- Keep suction bulbs of chest electrodes away from sunshine and excessive temperature;
- After long-term use, the surfaces of electrodes will be oxidized because of erosion and other causes. By this time, electrodes should be replaced to achieve high-quality ECG records.

8. Troubleshooting

8.1 ECG Waveform of Some Leads Cannot Be Displayed


- After the electrodes are connected to the patient, the ECG waveform is not stable and baseline drifts. The software is saturated or overflows when running.
- The equipment, patient cable and patient are not securely connected.
- The patient cable is faulty. Check the patient cable according to the method described in 4.3. If the problem still exists, please contact our after-sales service department or the designated repair unit.
- The signal channel of the equipment has malfunction. Please contact our after-sales service department or the designated repair unit.

8.2 Keypad Not Work

- Generally, it is resulted from the loose connector between the control panel and the circuit board. Please ask the professional maintenance personnel to open the cover of the equipment and reconnect the connector. If the problem still exists, please contact our after-sales service department or designated repair unit.

8.3 AC Interference

- In the ECG recording process, there is a certain range of regular interference. There is obvious ECG baseline jitter, as shown in the following ECG waveform.

This phenomenon may be caused by the following reasons, please check and eliminate in turn:

- Is the equipment grounded reliably?
- Are the electrodes or patient cable connected correctly?
- Are electrodes and skin applied with gel?
- Is the metal bed reliably grounded?
- Is the patient touching a metal part of the bed or wall?
- Has anyone else touched the patient?
- Are there any high-power electrical equipment nearby? (such as X-ray machines or ultrasonic instruments)
- Is the patient wearing jewelry such as glass or precious gemstone?
- The signal channel of the equipment has malfunction. Please contact our after-sales service department or the designated repair unit.

❗ If the AC interference cannot be removed after the above measures, please change the setting of AC Filter.

8.4 EMG Interference

The baseline jitter of the recorded ECG waveform is as follows :

This phenomenon may be caused by the following reasons, please check and eliminate in turn:

- Check whether the indoor environment is comfortable?
- Is the patient nervous?
- Is the bed too small for the patient?
- Did you talk to the patient during the recording?
- Are the electrode clamps on the limbs too tight?
- Are electrodes and skin applied with gel?
- The patient cable is faulty. Check the patient cable according to the method described in 4.3. If the problem still exists, please contact our after-sales service department or the designated repair point unit.
- The signal channel of the equipment has malfunction. Please contact our after-sales service department or the designated repair unit.

8.5 Baseline Drift

The baseline of the recorded ECG wave moves up and down irregularly, as shown below :

This phenomenon may be caused by the following reasons, please check and eliminate in turn:

- Is the electrode installation stable?
- Is the lead wire in good contact with the electrode?
- Are the electrodes and patient skin clean?
- Is gel applied to the skin where electrodes placed?
- Is it caused by movement or breathing of the patient?
- Are the old electrodes mixed with the new ones?

- The patient cable is faulty. Check the patient cable according to the method described in 4.3. If the problem still exists, please contact our after-sales service department or the designated repair unit.
- The signal channel of the equipment has malfunction. Please contact our after-sales service department or the designated repair unit.
- If the above measures still not working, please change the Baseline setting.

❗ If the above measures still not working, please change the Baseline setting.

8.6 Save ECG Data without Printing

This function can be enabled as long as the print paper is not installed in the equipment.

8.7 Paper-jam

- If it happens for the first time, the reason may be that the paper is not placed properly. Open the recording casing, pull out the paper from the carton, tear off the crumpled part, and then put it back into the casing. Close the casing after carefully adjusting the paper position.
- There is a problem with the print paper, replace it with a new one.
- If it is not the above situation, it may be the print module failure, please contact the manufacturer or local dealer for handling.

9. Warranty and After-sales Service

9.1 Warranty

The warranty scope: Failures caused by defects in materials and workmanship can be repaired or replaced for free during the warranty period.

- Main unit

Spengler warrants that ECGMAC's products meet the labeled specifications of the products. Under normal use and maintenance, the manufacturer can guarantee the free maintenance after receiving the report proving that there is a failure within one year from the date of delivery.

- Accessories :

Under normal use and maintenance, the manufacturer can guarantee the free maintenance of the accessories after receiving the report proving that there is a failure within six months from the date of delivery.

- Software :

The manufacturer can guarantee the free maintenance of the software after receiving the report proving that there is a failure within one year from the date of delivery. For software upgrades, please consult the manufacturer directly.

Note: Spengler's obligations under this warranty do not include freight and other expenses. All costs of repairing products beyond the warranty period shall be borne by the user. All maintenance and treatment shall be carried out by engineer and technical personnel approved by ECGMAC.

9.2 Manufacturer's Disclaimer

The warranty is void in cases of:

- Users disassemble and re-adjust the equipment by themselves
- Damage caused by alteration or repair by anyone not authorized by ECGMAC
- Replacement or removal of serial number label and manufacture label
- Damage caused by improper use
- Damage caused by abnormal use beyond the specified conditions

9.3 After-sale Service

If you have any question about maintenance, technical specifications or malfunctions of devices, please contact:

- the distributor: SPENGLER SAS - 30 rue Jean de Guiramand - 13290 Aix-en-Provence - FRANCE

- or manufacturer: Shenzhen ECGMAC Medical Electronics Co., Ltd

Address: 2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District, 518102 Shenzhen, China

Tel : 0755-27697821 ; 0755-27697823

Fax : 0755-27697823

Email : info@ecgmac.com

9.4 Service Life and Production Date

Production date: see product label for details

Service life of the product: 5 years (the service life is only limited to the main unit excluding the replaceable accessories; if the product needs to continue to be used after more than 5 years, it needs to be maintained every six months.)

Appendix A: Package and Accessories

A.1 Accessories

When the product is dispatched from the factory, the intact package shall contain the following accessories:

ECG Main Unit	1
Thermal Rolled Paper	1
Patient Cable	1
Chest Electrodes (suction ball)	6
Limb Electrodes (clamp)	4
Power Cord	1
Power Adapter	1
Grounding Wire	1
User Manual	1
User Manual of ECG Data Management Software (optional)	1
V1 CD of ECG Data Management Software (optional)	1

A.2 Caution

- 1) Open the box from the top;
- 2) After opening the box, check the accessories and user manual, and then check the instrument;
- 3) If the instrument does not work properly, contact the sales department or customer service department;
- 4) Please use the accessories supplied by ECGMAC. Accessories from other suppliers may damage the instrument and affect its performance and safety. Before using the accessories from other suppliers, please consult our customer service first;
- 5) To enable us to serve you in time, please fill out the warranty card (copy) and mail it to us;
- 6) The packing box should be kept properly for the use of regular verification or maintenance of the instrument.

Appendix B: Technical Specifications

B.1 Technical Specifications

1. Performance Specifications

Calibration Voltage	1mV $\pm$ 2%
A/D converter (resolution)	24bit (EM-301A/ EM-601A) 12bit (EM-301/ EM-301B/ EM-601/ EM-601B)
Input circuit current	$\leq 0.01 \mu A$
Sensitivity	$\pm 2\%$
Noise level	$\leq 15 \mu V$ (EM-301/ EM-301B/ EM-601/ EM-601B) $\leq 10 \mu V$ (EM-301A/ EM-601A)
DC offset voltage	$\pm 400 mV$ (EM-301/ EM-301B/ EM-601/ EM-601B) $\pm 1180 mV$ (EM-301A/ EM-601A)
Input impedance	$\geq 50M \Omega$ (10Hz)
Input voltage range	EM-301/ EM-301B/ EM-601/ EM-601B: each lead should no less than $(-7.5 \sim 7.5) mV$ EM-301A/ EM-601A: each lead should no less than $(-22.5 \sim 22.5) mV$
CMRR	EM-301/ EM-301B/ EM-601/ EM-601B: $\geq 100 dB$ (AC Filter off) ; $\geq 120dB$ (AC Filter on) EM-301A/ EM-601A: $\geq 110 dB$ (AC Filter off) ; $\geq 120 dB$ (AC Filter on)
Frequency response	0.05Hz $\sim$ 150 Hz (EM-301/ EM-301B/ EM-601/ EM-601B) 0.01Hz $\sim$ 350 Hz (EM-301A/ EM-601A)
Time constant	$\geq 3.2 s$ (EM-301/ EM-301B/ EM-601/ EM-601B) $\geq 5 s$ (EM-301A/ EM-601A)
Paper Speed	5 mm/s、6.25 mm/s、10 mm/s、12.5 mm/s、25 mm/s、50 mm/s ($\pm 2\%$)
HR range	30bpm $\sim$ 300bpm(± 1)
Gain	1.25 mm/mV、2.5 mm/mV、5 mm/mV、10 mm/mV、20 mm/mV、40 mm/mV、10/5 mm/mV、20/10 mm/mV、AGC. Tolerance $\pm 2\%$
Sampling Rate	2000 sps/channel (EM-301/ EM-301B/ EM-601/ EM-601B) 32000sps/channel (EM-301A/ EM-601A)
Storage rate	1000sps/channel
Pacemaker pulse detection	Amplitude: $\pm 2mV \sim \pm 700mV$ duration time: 0.1 ms $\sim$ 2.0 ms
Amplitude quantisation	4.563 μV /LSB (EM-301/ EM-301B/ EM-601/ EM-601B) 2.289 μV /LSB (EM-301A/ EM-601A)

2. Function Features

ECG input channel	Simultaneous 12 leads ECG acquisition
Background grid	Background grid display on the screen
Display color	Various color options of main interface background and waveform display
Lead off Alarm	When patient cable disconnected, there would be “lead off” alarm on display.
ECG waveform freeze and replay	Supports ECG waveform freeze and playback up to 300seconds
Support pacemaker detection	Support pacemaker detection, pacemaker detection sensitivity can be adjusted to low and high levels.
Working mode	Auto mode, Rhythm mode, Manual mode
Record mode	EM-601/EM-601A/EM-601B :3 4, 3 4+1R, 3 4+3R, 6 2, 6 2+1R, 1R EM-301/EM-301A/EM-301B : 3 4, 3 4+1R, 1R
Display mode	Simultaneous or Sequential
Sample mode	Pre-sample, real-time sample, trigger sample, Periodic sample
Date mode	ECG Mode, Demo Mode, Calibration Mode
Input mode	Floating ground input, with anti defibrillation protection circuit and pulse suppression function.
Lead system	Support Standard lead system, Posterior wall lead system, R-Pectoral lead system, R-pectoral/PW lead system, Up intercostal lead system, Down intercostal lead system, Cabrera lead system and User-defined lead name.
Defibrillation-proof Protection	defibrillation-proof Electric Shock Protection
Barcode Scanner Function	The patient information can be input by scanning one-dimensional and two-dimensional barcode.
System check function	It can self-check the machine in terms of display, keyboard, file system, battery, amplify board, network, and printer.
Patient records management function	Users can add or edit patient information. Equipped with file management function of patient records, including delete, upload, preview, import, export, query, print, etc.
Peripheral	Built-in SD card slot and USB port. Support SD card, U-disk, Keyboard, printer and barcode scanner
Local memory	Built in memory, save up to 1000ECG report, the data can be checked on PC.
Data Transmission	Equipped with wireless networking function, supports the transmission of ECG data to Bluetooth ECG and PC-ECG or ECG data management software, and supports the acquisition of appointment patient information from the server. (optional).
Report output format	Support PDF, JPEG, BMP, FDA-XML, EM-XML, SCP, BKG, DICOM and other formats output.
Print	1. Built-in Thermal Printer module. 2. It can be connected to designated external USB printer to print ECG report on A4 paper. 3. It supports print ECG grid on paper without grids.

3. Recorder

Recorder : Thermal Dot-matrix Recorder

Resolution	≥ 8 dots/mm (vertical direction) ≥ 40 dots/mm (25mm/s); ≥ 20 dots/mm (50 mm/s) (horizontal direction)
------------	---

Recorder Paper

Roll Paper : EM-601\EM-601A\EM-601B:110mm width
EM-301\EM-301A\EM-301B:80mm width

4. Display

7 inch TFT color display with 800 480 resolution

Background with grid, 12 leads ECG can be displayed on the same screen

Information can be displayed on LCD : No Paper, Lead Off Alarm, Operation Menu, Patient Information, ECG waveform, Battery Power, Date and Time, Heart Rate, 12 leads waveform, Working Mode, Paper Speed, Gain, and Filter.

5. Device Classification

1) Anti-electric-shock type: security category: EN 60601-1:2006, Class I & internally powered device.

2) Defibrillation-proof Type CF Applied part.

3) Degree of protection against harmful influx: Normal equipment (no protection against ingress of hazardous liquid).

4) Degree of safety of application in the presence of flammable gas: Not suitable for use under the presence of flammable gas.

5) Mode of operation: Continuous.

6) EMC: Group 1, Class A.

6. Others

12-lead patient cable with defibrillation proof.

Power supply

AC: 100 ~ 240 V, 50 / 60Hz, 0.17 A ~ 0.4 A

DC: 14.4 V / 5200 mAh, rechargeable battery

Power adapter output: DC19V 2.1A

Printer recommended: HP LaserJet 1020 (HP Company, US)

B.2 Physical Specifications

Main unit Dimensions : 310 mm 244 mm 65 mm

Package Dimensions : 380 mm 350 mm 190 mm

Net weight : 2.0 kg

Gross weight : 4.5 kg

B.3 Environment Condition

1. Transportation

Temperature	-20°C~+55°C
Relative Humidity	≤93%
Atmospheric Pressure	50kPa~106kPa

2. Storage

Temperature	-20°C~+55°C
Relative Humidity	≤93%
Atmospheric Pressure	50kPa~106kPa

3. Working

Temperature	+5°C~+40°C
Relative Humidity	≤80%
Atmospheric Pressure	86kPa~106kPa

Appendix C: Key Components

	Component	Model	Specifications	Note
1	Power adapter	LXCP40-019210	19V/2.1A	
2	Isolation transformer	EE16	DIP10	
3	LCD	AT070TN92		
4	Patient cables	ECG-FD08X4		

Appendix D: List of Functional Accessories

	Accessory	Model	Note
1	Patient Cable	ECG-FD08X4	
2	Limb Electrodes	ECG-FJX41	Optional
3	Chest Electrodes	ECG-FQX41	Optional
4	Limb Electrodes	ZJ-01	Optional
5	Chest Electrodes	XQ-01	Optional
6	Limb Electrodes	ZJ-02	Optional
7	Chest Electrodes	XQ-02	Optional

Appendix E: EMC Information

Instructions for use

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.

Warning : Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning : Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning : Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

Warning : Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment (detail model name), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If any : a list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g. shielded cable, load impedance) or specifically (e.g. by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).

If any : the performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term “ESSENTIAL PERFORMANCE” need not be used).

Technical description

- 1.all necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
2. Guidance and manufacturer’s declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacturer’s declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Applied

Table 2

Guidance and manufacturer’s declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	38 kV contact 32 kV, 34 kV, 38 kV, 315 kV air	38 kV contact 32 kV, 34 kV, 38 kV, 315 kV air
Electrical fast transient/ burst IEC 61000-4-4	Power supply lines : 32 kV input/output lines : 31 kV	Power supply lines : 32 kV
Surge IEC 61000-4-5	line(s) to line(s) : 31 kV. line(s) to earth : 32 kV. 100 kHz repetition frequency	line(s) to line(s) : 31 kV. 100 kHz repetition frequency
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 ° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	0% 0.5 cycle At 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 ° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz : 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150KHz to 80MHz : 3Vrms 6Vrms (in ISM bands) 80% Am at 1kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity							
	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27
	450	380 - 390	GMRS 460, FRS 460	FM 3.5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5240						
	5785						

Appendix F: Manufacturer Information

Manufacturer: Shenzhen ECGMAC Medical Electronics Co., Ltd

Factory Address: 2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District, 518102 Shenzhen, China

Tel: +86 755-27697821/ 27697823/ 27948579

Fax: +86 755-27697823-616

Website: www.ecgmac.com

Email : info@ecgmac.com

Appendix G: European Authorized Representative Information

Company Name: Well Kang Ltd

Company Address: Entreprise Hub, NW Business Complex, 1 Beraghmore Rd, Derry, BT54 8SE, Northern Ireland, UK

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
The screenshot shows the 'Information setup' screen. At the top, it says 'Information setup' on the left and 'Operate : Select Switch' on the right. The main area has two text input fields: 'Hospital name:' followed by a white rectangular box, and 'Device name:' followed by a white rectangular box. At the bottom, there are two buttons: 'Confirm F1' on the left and 'Cancel F2' on the right. A small 'en' button is located to the right of the 'Device name' field.

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The screenshot shows the 'Admin configuration' screen. At the top, it says 'Admin configuration' on the left and 'Operate : Select Switch' on the right. The screen contains several configuration options with radio buttons: 'ID generated by:' with options 'Acc.', 'Time' (selected), and 'User input'; 'Upload:' with options 'Auto' and 'Manual'; 'Upload Method:' with options 'FTP', 'Shared Directory', and 'DICOM'; 'File Format:' with options 'FDA-XML', 'EM-XML', 'SCP', 'BKG', and 'DICOM'; 'Picture Type:' with options 'Don't save', 'PDF', 'JPG', and 'BMP'; 'Picture Direction:' with options 'Horizontal' and 'Vertical'; 'Picture ECG Grid:' with options 'Off', 'Dot', and 'Line'; 'Picture Grid Color:' with options 'Gray', 'Blue', and 'Red'; 'Patient Query Protocol:' with options 'Off', 'HL7', and 'Worklist'; and 'Data Mode:' with options 'ECG mode', 'Demo', and 'Calibration'. At the bottom, there are five buttons: 'Language F1', 'Logo F2', 'Update F3', 'Barcode F4', and 'More F5'.

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File Format: ☐ FDA-XML ☐ EM-XML ☐ SCP ☐ BKG ☐ DICOM

Picture Type: ☐ Don't save ☐ PDF ☐ JPG ☐ BMP

Picture Direction: ☐ Horizontal ☒ Vertical

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