Spengler

LIAN® DESSILLONS

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1. PRODUCT IDENTIFICATION	12
2. INTENDED USE	14
3. CONTRAINDICATIONS, COMPLICATIONS AND CAUTIONS	15
4. TECHNICAL SPECIFICATIONS	15
5. MAINTENANCE AND CALIBRATION	16
6. USE AND PERFORMANCES OF THE DEVICE	17
7. CLEANING, STORAGE AND DISPOSAL	18
8. INFORMATION IN CASE OF SERIOUS INCIDENT	19
9. WARRANTY	19
10. TABLE OF USED SYMBOLS	19

The sphygmomanometer is a medical device that measures the pressure exercised by blood on arteries by the auscultatory method requiring a cuff or lower limb cuff and a stethoscope.

The sphygmomanometer is a mano-bulb type with aneroid diaphragm, very light.


The precision pressure gauge graduated from 0 to 300mmHg is fixed on an anodized aluminium valve, equipped with a decompression knob.

The readings of measured valued is facilitated by the size of the dial. Guaranteed pressure accuracy of $\pm 3\text{mm}$ over the whole measuring range.

- The medical class of the device is Im
- Lifetime is subjected to device use and cleaning and decontamination cycles. Usual lifetime is 2 years for the sphygmomanometer.
- This device is intended to be used with other devices such as: cuffs, lower limb cuffs and stethoscopes.
- The size of the cuff has to be adapted to patient's morphology.


LIAN® DESSILLONS references

DESIGNATION	REFERENCE	WITH SOFT CUFF®	WITH CLINIC® CUFF
Sphygmomanometer Lian® Dessillons Soft Cuff® small adult SPENGLER		TS50003S	
Sphygmomanometer Lian® Dessillons Soft Cuff® adult SPENGLER		TS50004S	
Sphygmomanometer Lian® Dessillons Soft Cuff® large adult SPENGLER		TS50005S	
Sphygmomanometer Lian® Dessillons Soft Cuff® multi cuffs SPENGLER		TS53000S	
Sphygmomanometer Lian® Dessillons Clinic® small adult cuff SPENGLER			TI50003S
Sphygmomanometer Lian® Dessillons Clinic® adult cuff SPENGLER			TI50004S
Sphygmomanometer Lian® Dessillons Clinic® large adult cuff SPENGLER			TI50005S
Sphygmomanometer Lian® Dessillons Clinic® multi cuffs SPENGLER			TI53000S
Sphygmomanometer Lian® Dessillons single manometer SPENGLER	A10370S		


References and sizes of the Soft Cuff® cuffs and lower limb cuff

REFERENCE WITH CONNECTOR	DESIGNATION	OUTPUT	CIRCUMFERENCE (CM)		COLOR TUBE
			MIN	MAX.	
BT101S	Newborn size bayonet cuff	1	6	12	BLACK
BT102S	Child size bayonet cuff	1	11	22	
BT103S	Little adult size bayonet cuff	1	16	28	
BT104S	Adult size bayonet cuff	1	26	33	
BT105S	Large adult size bayonet cuff	1	33	41	
CT104S	Lower limb adult size bayonet cuff	1	39	55	

References and sizes of the Clinic® cuffs and lower limb cuff

REFERENCE WITH CONNECTOR AND CAP	DESIGNATION	OUTPUT	CIRCUMFERENCE (CM)		COLOR TUBE
			MIN	MAX.	
BM101S	Newborn size bayonet cuff	1	6	12	BLACK
BM102S	Child size bayonet cuff	1	11	22	
BM103S	Little adult size bayonet cuff	1	16	28	
BM104S	Adult size bayonet cuff	1	26	33	
BM105S	Large adult size bayonet cuff	1	33	41	
CM104S	Lower limb adult size bayonet cuff	1	39	55	

INTENDED USE	Sphygmomanometers are designed for measuring non invasive blood pressure. The sphygmomanometer allows to measure the pressures (systolic and diastolic) of blood circulating in arteries. The readings are expressed in mmHg (millimeters of mercury) and include 2 values: • Systolic pressure • Diastolic pressure The blood pressure measurement is done with an adapted cuff to patient's morphology and a stethoscope. It then allows a medical diagnosis.										
INDICATION	Medical diagnosis, post-surgical surveillance, pregnant woman, patients with chronic disease (hypertension, diabetes, etc...). The non invasive measurement of blood pressure is done with a sphygmomanometer and a stethoscope. The sphygmomanometer allows to measure the back pressure on an artery coming from the inflated cuff while examining the artery.										
PATIENTS POPULATION	Newborns, pediatric patients and adults										
USERS	Users are healthcare professionals trained to the non invasive measurement of blood pressure technique using the auscultatory method. The device is not intended for use by profane people. User profil: <table data-bbox="314 877 955 1180"> <tr> <td>Occupation</td><td>Nurses, doctors, surgeons, mid-wives and other health care professionals.</td></tr> <tr> <td>Sex</td><td>Male/female</td></tr> <tr> <td>Age</td><td>>18</td></tr> <tr> <td>Level of education</td><td>Medical studies</td></tr> <tr> <td>Specific skills</td><td>Knowledge of blood pressure measurement technique</td></tr> </table>	Occupation	Nurses, doctors, surgeons, mid-wives and other health care professionals.	Sex	Male/female	Age	>18	Level of education	Medical studies	Specific skills	Knowledge of blood pressure measurement technique
Occupation	Nurses, doctors, surgeons, mid-wives and other health care professionals.										
Sex	Male/female										
Age	>18										
Level of education	Medical studies										
Specific skills	Knowledge of blood pressure measurement technique										

3. CONTRAINDICATIONS, COMPLICATIONS AND CAUTIONS

CONTRAINDICATIONS	No known contraindication
COMPLICATIONS	Side effect linked to the blood pressure measurement: none
CAUTIONS	<u>Cautions before use:</u> 1. Select the cuff adapted to patient's morphology 2. It is an absolute must that the cuff features the manometer adapted connector. 3. A check of the material is necessary before use: visually control the overall condition, and particularly the presence of the manometer, the decompression valve, the inflation bulb. If the device shows a dysfunction (cuff or manometer leakage), it has not to be used; <u>Cautions while using and after use:</u> 4. In case of device leakage, inflate again the cuff or change the whole device. 5. Cuff and lower limb cuff should not be manipulated with any sharp objects 6. Comply with cleaning and disinfection instructions, check the device after disinfection.

4. TECHNICAL SPECIFICATIONS

	MATERIAL	PRESSURE READING	FASTENING SYSTEM
Manometer	Aluminium body connector PVC bulb Polycarbonate glass Manometer protection: TPU Decompression valve: nickel-plated brass	Graduated manometer 0 - 300 mmHg, Precision ± 3 mmHg	Female (CMP001137) on manometer
SOFT CUFF® Cuff	PU leather Bladder in welded PU	N/A	Tube with A11313 (cuff male connector)
CLINIC® Cuff	PU coated fabric, high frequency welded	N/A	Tube with A11312 (cuff male connector)

- It is recommended to calibrate the manometer every 12 months.
- No liquid should enter the interior of the pneumatique cuff our pneumatic lower limb cuff, as it may partially or completely deteriorate the measurement device to which the cuffs are connected.
- In order to avoid risks, it is advised to use a DDM sealing system : the A11313 cap.
- Tubing should not be clamped with a clamp or an inappropriate device;
- Cuff and lower limb cuff should not be manipulated with any sharp objects.
- Cuff and lower limb cuff should not be inflated without being wrapped and closed with 2 velcro parts on a limb or on a sharpless rounded tubing .

List of spare parts: they have to be replaced by a biomedical technician.

Sphygmomanometer LIAN®DESSILLONS aluminium and polymer (PC)

REFERENCE	DESCRIPTION	SPECIFICATIONS
Soft Cuff® Cuff	Velcro cuff and bladder	Coated fabric and PU coated fabric for bladder
Clinic® Cuff	Monobloc cuff	PU coated fabric High frequency welded
A10362	Glass	PC
A10363	Upper manometer protector	TPU
A10369	Manometer housing	ABS PC
A10365	Body	Aluminium
A10226	Inflation bulb	PVC phthalate free
A10366	Lower manometer protection	TPU
A10367	0-300mmHg mechanism	Brass
A10368	Valve + filter + washer + circlips	
CMP001137	Female connector to screw	Nickel plated brass

INSTRUCTIONS FOR USE

1. It is an absolute must that the cuff features the manometer adapted connector.
2. Select the cuff adapted to patient's morphology Size and usage range are indicated on the cuff.
3. Leave the patient rest for 5 minutes prior to test.
4. Place the patient in a comfortable position, back and arm supported, feet uncrossed and flat on the floor without moving nor speaking.
5. Position the center of the bladder represented by a heart on the axis of the humeral artery.
6. Wrap the cuff properly around the arm and ensure its retention using the intended fastening system. The cuff should be positioned on the upper arm at the level of the right atrium.
7. Once the cuff is positioned and ready for use, the INDEX must be between the Min and the Max (heart axis). If the INDEX exceeds the Max index at the heart axis, replace the cuff with a larger size model. If the INDEX exceeds the Min index at the tip of the range, replace the cuff with a smaller size model.
8. When the cuff and INDEX are correctly positioned the blood pressure measurement process can begin.
9. Position the tips of the stethoscope lyre to your ears, the chestpiece in the listening position, just below the cuff, on the path of the artery.
10. Ensure that the pear pressure relief screw is in the closed position and inflate to 30 or 40 mmHg above the assumed pressure. The blood no longer circulates at the microphone, no sound is audible at the stethoscope.
11. Slightly unscrew the bulb decompression screw to slowly reduce the inflation pressure and allow blood to flow again, the first beats are noticeable from the stethoscope.
12. Read at this moment the pressure indicated on the manometer, this value corresponds to the maximum pressure or systolic pressure.
13. Continue decompression, the beats perceived with the stethoscope amplify and then fade to become inaudible.
14. Read at this time the indicated pressure, this reading corresponds to the minimum pressure or diastolic pressure.
15. Unscrew the decompression wheel further to completely purge the remaining air in the bladder.
16. Remove the cuff.

ESSENTIEL PERFORMANCE

The sphygmomanometer allows to measure the non invasive blood pressure with a cuff and a stethoscope

Cleaning

→ Manometer cleaning

- | | |
|----------|---|
| A | Spraying of decontaminant solution.
Do not use detergent, do not immerse |
|----------|---|

→ SOFT CUFF® & CLINIC® cleaning:

Recommended product : DETERQUAT / ANIOS.

- | | |
|----------|---|
| A | Spray the cuff with a detergent, disinfectant and bactericide foam, taking care to spread the product properly. |
| B | Leave on for 15 minutes |
| C | Rinsing is useless, wipe with a single use tissue |

→ SOFT CUFF® immersion cleaning:

Recommended product : ANIOS'CLEAN EXCEL D / STERANIOS

- | | |
|----------|--|
| A | Take off the bladder from the cuff. |
| B | Clean the cuff with water or in the washing machine at 40° (60 ° being the maximum), without spinning. The number of machine washes as well as the washing temperature reduces the duration of use of the cuff |

→ CLINIC® immersion cleaning:

Recommended product : ANIOS'CLEAN EXCEL D / STERANIOS

- | | |
|----------|---|
| A | Seal the tube(s) with the closure system provided: - Connector with cap: A11313 |
| B | Immerse the tube and cuffs in a solution. |
| C | Immersion time: 15 minutes. |
| D | Rinse thoroughly with water, renewal of the soaking bath after each use. |

Secured disposal of the device

The device should be disposed off, after being decontaminated in compliance with current local practices.

It should be discarded in a specific disposal container.

Usage conditions

TEMPERATURE	HYGROMETRY	PRESSURE
+ 10°C to + 40°C	15% à 85% (without condensation)	700 hPa to 1060 hPa

Storage conditions

PACKAGING	TEMPERATURE	HYGROMETRY	PRESSURE
Origin	-20°C à 70°C	15% à 90% (without condensation)	700 hPa à 1060 hPa

Any serious incident with the device should be reported and notified to the manufacturer and to the competent Health Authority of the State Member where the patient lives.

9. WARRANTY

WARRANTY	Spengler replaces or repairs any non-functioning device as notified in the manual. The sphygmomanometer is guaranteed for 24 months. The cuffs are guaranteed for 12 months.
WARRANTY LIMIT	the product has to be used in compliance with the instructions for use and with the foreseen indications . It should not be modified neither accidentally deteriorated before use.

10. SYMBOLES USED

SYMBOL	INFORMATION	SYMBOL	INFORMATION
	Manufacturer		Temperature limits
	Medical Device		Humidity limits
	Reference		Atmospheric pressure limits
	Manometer		Pouch
	Manufacturing date		CE
	Lot number		Consult instructions for use
	Serial number		Swiss representative
	Cuff		User manual

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DDM appartient à Spengler SAS
Spengler S.A.S | 30 rue Jean de Guiramand |